## Instructions for Use

# BioMend® and BioMend® Extend™ Absorbable Collagen Membrane

**Distribuido by:**  
Zimmer Biomet Dental  
4555 Riverside Drive  
Palm Beach Gardens  
FL 33410 USA  
+1.561.776.6700  
zimmerbiometdental.com

## ENGLISH

### Instructions for Use

## BioMend® and BioMend® Extend™ Absorbable Collagen Membrane

Before using this product, the surgeon/practitioner should carefully study the indications, contraindications, recommendations, warnings and instructions and fully comply with them. The manufacturer, the importer and the distributor of these products are not liable for complications, other negative effects or damages that might occur for reasons such as incorrect indications or surgical technique, unsuitable choice of material or handling thereof, unsuitable use or handling of the instruments, asepsis and so on. The surgeon/practitioner is responsible for any such complications or other consequences. It is also the surgeon/practitioner's responsibility to properly instruct and inform the patient on the functions, handling and necessary care of the product and on all known product risks.

**Description**  
*BioMend* and *BioMend Extend* membranes are white, compressed, non-friable matrices fabricated from collagen derived from bovine deep flexor (Achilles) tendon. Bovine tendon is known to be one of the purest sources of Type I collagen that can be readily obtained and processed in commercial amounts. *BioMend* and *BioMend Extend* membranes are completely absorbable, eliminating the need for the second surgical procedure often required to remove a non-resorbable membrane.

Under scanning electron microscopy, *BioMend* and *BioMend Extend* absorbable collagen membranes have a morphology of dense oriented fibers for mechanical strength. Macromolecular permeation studies have shown that the membrane is permeable to macromolecules. *BioMend* and *BioMend Extend* membranes can be cut to any size or shape in the wet or dry state without tearing or fragmenting.

*BioMend* and *BioMend Extend* membranes have an effective pore size of 0.004 microns, which is designed to help retard epithelial down growth during early phases of healing. Being semi-occlusive, the membrane allows essential nutrients to pass through it. *BioMend* and *BioMend Extend* membranes incorporate into the surrounding tissue. The *BioMend* membrane should be absorbed approximately 8 weeks following surgery. The *BioMend Extend* membrane should be absorbed approximately 18 weeks following surgery.

*BioMend* and *BioMend Extend* membranes are sterilized using ethylene oxide gas.

**Indications**  
*BioMend* and *BioMend Extend* absorbable collagen membranes are absorbable, implantable materials indicated for:

- Guided tissue regeneration procedures in periodontal defects to aid in wound healing post periodontal surgery.
- Use in dental surgery procedures as a material for placement in the area of an implant, bone defect or ridge construction to aid in wound healing post dental surgery.

**Contraindications**

- BioMend* and *BioMend Extend* absorbable collagen membranes are contraindicated in patients who have acute infections or contaminated wounds in the oral cavity.
- BioMend* and *BioMend Extend* absorbable collagen membranes are contraindicated in patients with a known history of allergic responses to collagen.
- BioMend* and *BioMend Extend* absorbable collagen membranes are contraindicated in patients who are allergic to bovine-derived products.
- BioMend* and *BioMend Extend* absorbable collagen membranes are contraindicated in patients who have clinically significant renal, hepatic, cardiac, endocrine, hematologic, autoimmune or systemic disease, which, in the physician's judgment, will prevent safe implantation or likely healing.

**Warnings**  
Clinicians should use care in screening their patients for any known allergies to collagen or bovine-derived products.

Hypersensitivity reactions have been noted with the use of other products containing bovine collagen; therefore, the possibility exists of developing a local sensitivity response to *BioMend* and *BioMend Extend* membrane.

**Precautions**  
It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects, and precautions as well as to seek the services of a trained dental professional if there are any changes in the performance of the membrane (e.g., infection, pain, any other unusual symptoms that the patient has not been told to expect). If these conditions occur, the patient should be instructed to see a trained dental professional immediately.

As with all surgical procedures, caution should be exercised when treating medically compromised patients such as patients receiving long-term steroidal therapy or currently taking anticoagulants.

Patients with clinically significant diseases, indicating a history of anaphylactic reactions, autoimmune disease, uncontrolled diabetes or severe hypertension have not been implanted with the device; therefore the safety and effectiveness for these patients have not been demonstrated.

The safety and effectiveness of the device has not been evaluated in pregnant women or children. Therefore, caution should be used in these patients.

*BioMend* and *BioMend Extend* membranes cannot be resterilized. Opened, unused *BioMend* and *BioMend Extend* membrane must be discarded. In vivo stability may be adversely affected if resterilized. Cross contamination and infection may occur if reused.

Do not use if the product sterilization barrier or its packaging is compromised.

*BioMend* and *BioMend Extend* absorbable collagen membranes are not intended for use on defects outside the indications stated.

*BioMend* and *BioMend Extend* absorbable collagen membranes have not been clinically evaluated in patients with conditions involving extremely severe defects with little remaining periodontium.

**Adverse Reactions**  
Possible complications which can occur with any dental surgery include swelling of the intraoral tissue, thermal sensation, gingival recession, gingival bleeding, flap sloughing, resorption or ankylosis of the treated root, some loss of crestal bone height, infection, pain or complications associated with the use of anesthesia.

As with any type of surgical therapy, the patient may experience minor discomfort for a few days.

Spontaneous exfoliation of the material may occur in the immediate postoperative period if *BioMend* and *BioMend Extend* membranes are not adequately covered by the mucogingival flap.

**Sterility**  
*BioMend* and *BioMend Extend* membranes are supplied sterile, non-pyrogenic and for single use only. The absorbable collagen membrane should be removed from their packaging using sterile gloves or instruments.

**Administration**  
It is advisable that only clinicians trained in related treatment and planning and in the technique of placing dental membranes should use *BioMend* and *BioMend Extend* absorbable collagen membranes.

Mucoperiosteal incision flaps are developed in the site to be treated. The incision should be sulcular when possible. The clinician should perform thorough debridement. As much tissue as possible should be preserved to allow for primary closure of the wound and correct positioning of the flaps. Space-making material such as autologous bone, mineralized bone particles, demineralized bone matrix and inorganic grafting materials may be used to fill the defect.

*BioMend* and *BioMend Extend* membranes can be placed either dry or hydrated. If the clinician prefers the handling characteristics of the hydrated collagen, the membrane can be hydrated in sterile water or saline solution prior to final placement.

The product can be trimmed to the size and shape of the defect in the dry or wet state using sharp, sterile scissors.

*BioMend* and *BioMend Extend* membranes are placed over the defect, in oral surgical procedures, and as close to the tooth as possible. The absorbable collagen membrane should extend a minimum of 2-3 mm beyond the bony defect apically, mesially and distally. Additional trimming of the membrane may be performed, avoiding sharp corners that could perforate overlying tissue. *BioMend* and *BioMend Extend* membranes can be sutured in place if desired using absorbable sutures and a non-cutting needle. Resorbable tacks can also be used to affix the membrane. Gingival flaps should be coronally positioned over the defect and the absorbable collagen membrane. The mucoperiosteal flap should completely cover the absorbable collagen membrane if possible and be sutured in place.

**Postoperative Procedures**

*BioMend* and *BioMend Extend* absorbable collagen membranes are fully absorbable and should not be removed.

Periodontal packing may be applied to the wound site. While this is a matter of clinician preference, care must be taken not to over compress the area. Post-operative care should include the following minimum steps:

Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate twice daily for four weeks following surgery. The wound site may additionally be swabbed with a cotton-tipped applicator dipped in the antimicrobial agent.

The patient should refrain from brushing the treated area for two weeks following surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Instruction will be dependent on an evaluation of wound healing. Dental floss should not be used prior to four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated.

The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to six to eight weeks following surgery. The patient may then return to a normal oral hygiene routine.

The *BioMend* membrane should be absorbed approximately 8 weeks following surgery. The *BioMend Extend* membrane should be absorbed approximately 18 weeks following surgery. However, probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature regenerated tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

**Storage**  
The product should be stored at room temperature. Avoid excessive heat and humidity.

**How supplied**  
Individually packaged in a variety of sizes:

15 mm x 20 mm, 1/box; 20 mm x 30 mm, 1/box; 30 mm x 40 mm, 1/box

**Caution**

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician or dentist.

U.S. Patent #7,807,192

©2017 Zimmer Biomet. All rights reserved.

BioMend and BioMend Extend are registered trademarks of Zimmer Biomet or its affiliates.

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