# NOVABONE®

## Dental Putty

### Bioactive Synthetic Graft

#### Grefte synthétique bioactive


#### Bioaktives Synthetisches Transplantate

#### Injerto sintético bioactivo

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#### Enxerto sintético bioactivo

#### Bioaktief synthetisch implantaat

# NOVABONE®

## Dental Putty

### Bioactive Synthetic Graft

#### Instructions For Use

##### Indications for Use:

The intended use of NovaBone Dental Putty is to provide a safe, biocompatible, synthetic bone graft material for use in oral, dental intraosseous, and craniofacial defects. It is used alone in a manner comparable to autogenous bone graft chips or allograft bone particulate (Demineralized Freeze Dried Bone) or may be mixed with either (typically 1:1 ratio v/v) as a bone graft extender. Typical uses include:

- Periodontal/Infrabony defects,
- Ridge Augmentation (sinusotomy, osteotomy, cystectomy),
- Extraction sites (ridge maintenance/augmentation, implant preparation/placement),
- Sinus lifts
- Cystic cavities
- Craniofacial augmentation

For larger defects, a mixture of NovaBone Dental Putty with an equal volume of allograft or autograft bone and bone marrow may improve new bone formation.

##### Description:

NovaBone Dental Putty is an osteoconductive bioactive device used for alloplastic bone grafting of osseous defects. It is a pre-mixed composite of bioactive calcium-phospho-silicate particulate and a synthetic, absorbable binder. The bioactive particulate is composed solely of elements that exist naturally in normal bone (Ca, P, Na, Si, O). The absorbable binder is a combination of polyethylene glycol and glycerin. The device requires no mixing or preparation prior to application. The non-hardening putty is supplied ready-to-use, to be applied directly to the intended graft site. The binder is then absorbed from the site such that only the bioactive particulate remains. NovaBone Dental Putty is partially radiopaque, which enables the clinician to make immediate, as well as long-term, post-operative evaluation. Positive surgical results are related directly to factors of patient cooperation, diligent home care practices, proper case selection and surgical technique.

##### Contraindications:

NovaBone Dental putty should not be used in patients with:

1. Are taking immunosuppressant medications, or other medications known to affect the skeleton (e.g. chronic glucocorticoid usage >10mg/day for the previous 3 months). Estrogen replacement therapy is allowed.
2. Need chronic anticoagulant therapy (e.g. heparin). Prophylactic use of Coumadin or aspirin postoperatively is allowed.
3. Have a systemic metabolic disorder known to adversely affect mucosla healing or bone healing and mineralization (e.g. poorly controlled insulin-dependent diabetes, renal osteodystrophy, Paget's disease), other than primary osteoporosis.
4. Have had, or are undergoing irradiation treatment of the graft region.

In addition, the prognosis in periodontal applications must be considered less positive if home care is inadequate, if endodontic or pulpal problems exist, if the patient is on a steroid regimen that may cause bone destruction, or where mobile teeth are treated.

##### Instructions for Use:

Routine surgical procedures should be used to expose the surgical site. Once exposed, eliminate all granulation or necrotic tissue at the defect site. Irrigate the defect with sterile saline or water and evacuate the excess. NovaBone Dental Putty requires no special handling or mixing procedures prior to use. All device packaging should be inspected prior to use to insure maintenance of sterility.

1. Remove device from sterile packaging.
2. For cartridge system, snap the cartridge into place in the Dispenser Handle, and remove the black cap from the cannula.
3. Express the desired amount either into a sterile dish or into the defect site.
4. If autograft or allograft bone is to be used, ensure that it is fully hydrated prior to mixing with an amount approximately equal to the NovaBone Dental Putty to be used. Mixing time is not critical.
5. NovaBone Dental Putty will remain in a pliable condition until application. NovaBone Dental Putty does not set like a cement.
6. After placement of the NovaBone Dental Putty, remove any excess material and close as per standard practice. Ensure primary closure of the soft tissues over the graft site if required a surgical dressing may be placed over the defect for 1-2 weeks.
7. The device cannot be reused after implantation due to risk of device or patient contamination. Excess material cannot be saved for later use due to risk of contamination and infection; all excess material and packaging must be discarded.
8. Postoperative antibiotic, analgesic and home care regimen should be prescribed as required.

##### Warnings:

Possible complications are the same as to be expected of autogenous bone grafting procedures. Complications that may arise as a result of the surgical procedure may include: oral indications - tooth sensitivity, gingival recession, flap sloughing, resorption or ankylosis of the treated root, abscess formation; general indications – pain, swelling, superficial wound infection, deep wound infection, deep wound infection with osteomyelitis, loss of reduction, loss of bone graft, graft protrusion and / or dislodgement, and general complications associated with anesthesia use and / or surgery.

NovaBone Dental Putty does not possess sufficient mechanical strength to support load bearing defects prior to soft and hard tissue ingrowth. In cases of use in load bearing regions such as mandibular fractures, standard internal or external stabilization techniques must be followed to obtain rigid stabilization in all planes.

NovaBone Dental Putty is intended for manual application and is not intended for injection through a constrained opening or under high pressure. The syringe packaging must not be modified to permit high pressure delivery of the device. High pressure injection of NovaBone Dental Putty should not be conducted as it could result in device over-pressurization, which may lead to device extrusion beyond the intended application site or to embolization of fat or the device into the bloodstream.

##### Precaution:

NovaBone Dental Putty is not intended for use in defects outside the indication statement and has not been clinically tested for use in pregnant women. Existing pathological conditions (e.g. infections) may compromise the results and should be controlled or eliminated prior to the use of NovaBone Dental Putty.

##### Material Notes – Osteostimulation:

NovaBone Dental Putty is an osteoconductive bone graft device. *In vivo* tests have demonstrated more extensive bone formation at early post-implantation time periods than for other standard osteoconductive devices such as hydroxyapatite.<sup>1,2</sup> *In vitro* cell culture tests with both animal and human osteoblasts have demonstrated an osteostimulative effect, defined as the active stimulation of osteoblast proliferation and differentiation as evidenced by increased levels of DNA synthesis and of the osteoblast markers osteocalcin and alkaline phosphatase.<sup>3-6</sup> This stimulation has been attributed as being the result of the interaction between osteoblasts and the ionic dissolution products released from NovaBone Putty particles during their absorption. Clinical data on this acceleration of bone formation and the increased levels of DNA synthesis, osteocalcin and alkaline phosphatase have not been established in humans.

##### Stability:

The device is provided STERILE unless the package is open or damaged. Do not use if the sterile packaging is damaged. The content of each package is designed for single use only. Device cleaning and resterilization is not possible. Shelf life is 4 years. Do not use after expiration date. Store in a dry place (<25°C).

##### Caution:

Federal law limits the device to sale by, or on the order of, a licensed dentist or physician.

##### Presentation:

- Pre-mixed device in single-use package sizes:

	
<b>Manufacturer:</b> NovaBone Products, LLC 13510 NW US Highway 441 Alachua, Florida 32615 USA Tel: 1-386-462-7660 www.novabone.com	<b>EU Authorized Representative:</b> Peter Smith 18 Yeates Close Thame, Oxfordshire, OX9 3AR UK Tel: (44) 184 421 6243

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