## CONTENT

<b>EN (English)</b> Instruction Manual Zfx™ Hex / Hexalobular Screwdriver.....	2
<b>BG (Български)</b> Ръководство за употреба Шестостенна / шестоъчна отвертка Zfx™.....	5
<b>CS (Česky)</b> Návod k použití Šestihřanný / hexalobulární šroubovák Zfx™.....	8
<b>DA (DANISH)</b> Brugervejledning Zfx™ Hex / Hexalobular Screwdriver.....	11
<b>DE (Deutsch)</b> Gebrauchsanweisung Zfx™ Sechskant- / Sechsrund-Schraubendreher.....	14
<b>EL (Ελληνικά)</b> Εγχειρίδιο οδηγιών Zfx™ Εξάγωνο/εξάλοβο κατσαβίδι.....	17
<b>ES (Español)</b> Manual de instrucciones Destornillador Zfx™ Hex / Hexalobular.....	20
<b>ET (Eesti)</b> Kasutusjuhend Zfx™ Hex / Hexalobular kuuskant-kruvikeeraja.....	23
<b>FI (Suomi)</b> Käyttöopas Zfx™ kuusiokolo / kuusikulmainen ruuvimeisseli.....	26
<b>FR (Français)</b> Mode d'emploi Tournevis hexagonal / hexalobé Zfx™.....	29
<b>HR (Hrvatski)</b> Priručnik Zfx™ Šesterokutni / Heksalobularni Odvijač.....	32
<b>HU (Magyar)</b> Kezelési útmutató Zfx™ imbusz / torx csavarhúzó.....	35
<b>IT (Italiano)</b> Manuale di istruzioni Cacciavite esagonale / esalobato Zfx™.....	38
<b>LT (Lietuvių k)</b> Instrukcijų vadovas „Zfx™“ šešiakampis / šešiaskiltis atsuktuvus.....	41
<b>LV (Latviski)</b> Lietošanas rokasgrāmata Zfx™ sešstūrveida /heksalobulārais skrūvgriezis.....	44
<b>NL (Nederlands)</b> Instructiehandboek Inbus- / hexalobulaire Zfx™-schroevendraaier.....	47
<b>PL (Polski)</b> Instrukcja obsługi Zfx™ Wkręтак sześciokątny / torx.....	50
<b>PT (Português)</b> Manual de Instruções Chave de parafusos Hexagonal / Hexalobular Zfx™.....	53
<b>RO (Română)</b> Manual de instrucțiuni Șurubelniță hexagonală /hexalobată Zfx™.....	56
<b>SK (Slovak/Sloavkian)</b> Návod na použitie Zfx™ Skrutkovač na skrutky s vnútorným šesťhranom / hexalobulárom.....	59
<b>SL (Slovenščina)</b> Navodila za uporabo Šestrobi / heksalobularni izvijač Zfx™.....	62
<b>SV (Svenska)</b> Instruktionsmanual Zfx™ Hex / Hexalobulär Skruvmejsel.....	65

# Instruction Manual


## Zfx™ Hex / Hexalobular Screwdriver

Version: 1/06.2021  
Compiled at: 08.06.2021  
Effective date: 08.06.2021  
Replaces version: -

### Indications for Use

The Hex/Hexalobular Screwdriver is used to tighten/loosen Hex/Hexalobular (laboratory try-in or final prosthetic) Screws to enable assembly/disassembly of single-unit (engaging) and multi-unit (non-engaging) restorative systems containing the Hex/Hexalobular screw to corresponding Zimmer Biomet Dental Implant Systems.

HLD –Hexalobular Drive

HD –Hex Drive

### Device description

The Zfx™ Hex / Hexalobular Screwdriver is used to drive the Hex/Hexalobular screw in a straight or angulated position. The limitation of angulation for Hexalobular is 25° and Hexagonal 20°. The Screwdriver are available in different overall length per size, refer to table 1.

The screwdriver is designed with an ISO-latch connection allowing for use with Torque Wrench or appropriate other fabricates adapters according the ISO 1797 Standard.

The screwdriver is an individually packaged reusable component and is available in various lengths. The screwdriver is provided non-sterile and must be cleaned and sterilized prior to use in a patient's mouth.

### Materials

The Zfx™ Hex / Hexalobular Screwdriver are manufactured from medical grade stainless steel.


### Procedure

#### × Handling

1. Connect the ISO-latch portion of the driver into the Torque Wrench or the Hand Piece.
2. Engage the correct screwdriver to the corresponding Hexagonal / Hexlobular Screw to tighten or loosen.

Refer to the table on page 3 for maximum torque recommendations for the use of Zfx™ Hex / Hexalobular Screwdrivers.

**Please note:** For the recommended final Prosthetic Screw torque values please refer to the label on the screw package or to product catalog information!

#### × Cleaning and Sterilization

The Hex/Hexalobular Screwdriver is non-sterile. Prior to use, the driver must first be cleaned and sterilized. The driver can be used up to 15 cleaning/sterilization cycles.

#### Cleaning

1. Rinse the driver under running cold tap water for 2 minutes.
2. Hand wash instruments in a bath with demineralized water for 5 minutes. Remove debris with a soft brush where possible. Cavities should be rinsed with a syringe. This should be repeated until all debris is removed.
3. Clean the instrument in an ultrasonic bath using a mildly alkaline cleaner (1.5%) at a frequency of at least 40 kHz for 20 minutes.
4. Rinse thoroughly with demineralized water for 1 min.

#### Sterilization

The recommended sterilization procedure is full cycle pre-vacuum steam sterilization at a temperature of 134°C (273.2° F) for 3 minutes, 4 pulse minimum or full cycle gravity steam sterilization at 134°C (273.2° F) for 3 minutes minimum. The products are ready for use 5 minutes after completion of the autoclaving process.

**Note:** Follow the loading instructions issued by the sterilizer manufacturer.

#### Storage and Handling

The Hex/Hexalobular Screwdriver should be stored at room temperature.

#### Contraindications

None

#### Known Side Effects

None

# Instruction Manual

## Zfx™ Hex / Hexalobular Screwdriver

Version: 1/06.2021  
Compiled at: 08.06.2021  
Effective date: 08.06.2021  
Replaces version: -

### Warnings

- × The Zfx™ Hex/Hexalobular Screwdriver must only be used and handled by dental professionals who have approved training in dentistry.
- × The Zfx™ Hex / Hexalobular Screwdriver may only be used with the corresponding Screw size (refer to table 1).
- × The Zfx™ Hex / Hexalobular Screwdriver are small and, therefore, they must be handled with caution to avoid swallowing or inhaling by the patient.
- × The use of a screw torque higher than the recommended value may result in damage to the Screwdriver, TiBase, screws, custom abutment, and / or the implant.
- × The Zfx™ Hex / Hexalobular Screwdriver must be in good condition. A worn screwdriver can lead to over-rotation of the hexagonal / hexalobe in the screw head which can result in the screw not being fully tightened and / or not being able to be removed.

### Explanation of symbols

Non-sterile

Do not use if package is damaged

Manufacturer Information: Address of legal manufacturer facility

Unique Device Identifier

LOT-Number, Symbol followed by Production Date of the lot.

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

Quantity / Unit of package content

Consult instructions for use  
[www.zfx-dental.com](http://www.zfx-dental.com)

Consult instructions for use for mating Hex/Hexalobular screws and angulated screw channel TiBase abutments

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Refer to the below table for maximum torque recommendations for the use of Zfx™ Hex / Hexalobular Screwdrivers.

**Table 1**

Item Name	Article No.	Max. Torque (Ncm)	Driver Size
Zfx™ Hexalobular Screwdriver EIGER, 21 mm	ZFX02HLD15421	40	HLD 1.54
Zfx™ Hexalobular Screwdriver EIGER, 28 mm	ZFX02HLD15428	40	HLD 1.54
Zfx™ Hexalobular Screwdriver EIGER, 36 mm	ZFX02HLD15436	40	HLD 1.54
Zfx™ Hexalobular Screwdriver JORASSES, 21 mm	ZFX02HLD16821	40	HLD 1.68
Zfx™ Hexalobular Screwdriver JORASSES, 28 mm	ZFX02HLD16828	40	HLD 1.68
Zfx™ Hexalobular Screwdriver JORASSES, 36 mm	ZFX02HLD16836	40	HLD 1.68
Zfx™ Hexalobular Screwdriver MERU, 21 mm	ZFX02HLD18021	40	HLD 1.8
Zfx™ Hexalobular Screwdriver MERU, 28 mm	ZFX02HLD18028	40	HLD 1.8
Zfx™ Hexalobular Screwdriver MERU, 36 mm	ZFX02HLD18036	40	HLD 1.8
Zfx™ Hex Screwdriver AORAKI, 21 mm	ZFX02HD12021	35	HD 1.2
Zfx™ Hex Screwdriver AORAKI, 28 mm	ZFX02HD12028	35	HD 1.2
Zfx™ Hex Screwdriver AORAKI, 36 mm	ZFX02HD12036	35	HD 1.2
Zfx™ Hex Screwdriver BADILE, 21 mm	ZFX02HD13321	40	HD 1.33
Zfx™ Hex Screwdriver BADILE, 28 mm	ZFX02HD13328	40	HD 1.33
Zfx™ Hex Screwdriver BADILE, 36 mm	ZFX02HD13336	40	HD 1.33
Zfx™ Hex Screwdriver DENALI, 21 mm	ZFX02HD14021	40	HD 1.4
Zfx™ Hex Screwdriver DENALI, 28 mm	ZFX02HD14028	40	HD 1.4
Zfx™ Hex Screwdriver DENALI, 36 mm	ZFX02HD14036	40	HD 1.4
Zfx™ Hex Screwdriver FUJI, 21 mm	ZFX02HD12721	40	HD 1.27
Zfx™ Hex Screwdriver FUJI, 28 mm	ZFX02HD12728	40	HD 1.27
Zfx™ Hex Screwdriver FUJI, 36 mm	ZFX02HD12736	40	HD 1.27
Zfx™ Hex Screwdriver PISCO, 21 mm	ZFX02HD15021	40	HD 1.5
Zfx™ Hex Screwdriver PISCO, 28 mm	ZFX02HD15028	40	HD 1.5
Zfx™ Hex Screwdriver PISCO, 36 mm	ZFX02HD15036	40	HD 1.5
Zfx™ Hex Screwdriver VINSON, 21 mm	ZFX02HD16021	40	HD 1.6
Zfx™ Hex Screwdriver VINSON, 28 mm	ZFX02HD16028	40	HD 1.6
Zfx™ Hex Screwdriver VINSON, 36 mm	ZFX02HD16036	40	HD 1.6
Zfx™ Hex Screwdriver EVEREST, 28 mm	ZFX02HD10028	15	HD 1.0

4099EN REV/06/21 © 2021 Zfx GmbH & Zimmer Dental Inc. | Zimmer and the Zimmer logo are trademarks of Zimmer Inc. or its affiliates and Zfx and the Zfx logo are trademarks of Zfx GmbH, Germany. All other trademarks are the property of the respective owner.

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