## RegenerOss® Resorbable Xenograft Porcine Anorganic Bone Mineral

### Instructions for Use

**Manufacturer:** Collagen Matrix, Inc. 15 Thornton Road Oakland, New Jersey 07436 USA www.collagenmatrix.com

**Distributed By:** Zimmer Biomet Dental 4555 Riverside Drive Palm Beach Gardens, FL 33410 USA +1.561.776.6700 zimmerbiometdental.com

## ENGLISH

**RegenerOss® Resorbable Xenograft Porcine Anorganic Bone Mineral**

#### INSTRUCTIONS FOR USE

A sterile, biocompatible anorganic porous bone mineral for use in periodontal, oral and maxillofacial surgery.

#### Description:

*RegenerOss® Resorbable Xenograft* is a porous bone mineral matrix consisting predominantly of calcium phosphate. It is produced by removal of the organic components from porcine cancellous bone. *RegenerOss® Resorbable Xenograft* is sterilized by gamma irradiation. The product is available in granular form, non-pyrogenic and for single use only.

#### Properties/Actions:

The anorganic bone matrix of *RegenerOss® Resorbable Xenograft* has interconnecting macro- and microscopic porous structure that supports the formation and ingrowth of new bone at the implantation site. The use of *RegenerOss® Resorbable Xenograft* may be considered when autogenous bone is not indicated, or insufficient in quantity to fulfill the needs of the proposed surgical procedure.

#### Indications and Usage:

*RegenerOss® Resorbable Xenograft* is indicated for:

- Augmentation or reconstructive treatment of the alveolar ridge.
- Filling of infrabony periodontal defects.
- Filling of defects after root resection, apicoectomy, and cystectomy.
- Filling of extraction sockets to enhance preservation of the alveolar ridge.
- Elevation of the maxillary sinus floor.
- Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR).
- Filling of peri-implant defects in conjunction with products intended for Guided Bone Regeneration (GBR).

#### Instructions for Use:

- After exposure of the bony defect with a mucoperiosteal flap, all granulation tissue must be carefully removed.
- Mix *RegenerOss® Resorbable Xenograft* with autogenous bone, osseous coagulum, patient’s blood or sterile normal saline. If large maxillofacial defects are present, *RegenerOss® Resorbable Xenograft* should be mixed with autogenous bone in a ratio of approximately 1:1.
- In order to assure the formation of new bone, *RegenerOss® Resorbable Xenograft* should only be placed in direct contact with well vascularized bone. Cortical bone should be mechanically perforated.
- Loosely pack *RegenerOss® Resorbable Xenograft* into osseous defect using a sterile instrument. Use of excessive force will result in compression of the particles and loss of trabecular architecture.
- Overfilling of the defects should be avoided.
- The mucoperiosteal flaps should be sutured to achieve primary closure, if possible. A surgical dressing may be placed over the wound for one to two weeks.
- If primary wound closure cannot be achieved completely, further immobilization of the flap (e.g., by incision through the periosteum) should be performed and/or a bioabsorbable membrane (e.g. Resorbable Collagen Membrane) should be placed over the bone graft site.
- Post-operative procedures: Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for 2 weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent. Consult the manufacturer’s Instructions for Use for side effects and other known risks associated with the use of chlorhexidine gluconate oral rinses.

#### Contraindications:

Contraindications customary to the use of bone grafts should be observed. *RegenerOss® Resorbable Xenograft* should not be used in patients with:

- Acute or chronic infection (osteomyelitis) at the surgical site
- Metabolic diseases (diabetes, hyperparathyroidism, osteomalacia)
- Severe renal dysfunction, severe liver disease
- High dose corticosteroid therapy
- Vascular impairment at the implant site
- Osteoporosis
- Known allergy or hypersensitivity to porcine-derived implantable materials

#### Warnings:

The device should be secured to prevent motion and migration; use in areas where the graft can be adequately contained.

Do not use if package is opened or damaged or if expiration date has been exceeded.

*RegenerOss® Resorbable Xenograft* cannot be re-sterilized or re-used. Open, unused product must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

#### Precautions:

In order to facilitate the formation of new bone, *RegenerOss® Resorbable Xenograft* should only be implanted in direct contact with a well vascularized bony tissue. Drilling may be recommended to facilitate bleeding from cortical bone.

In larger defects a mixture of autogenous bone or bone marrow may improve the formation of new bone.

#### Implantology

Generally, in augmented areas, the placement of titanium fixtures should take place once the bone has sufficient strength and integrity for dental implant placement, which is typically greater than 6 months after implantation of a bone graft material. For sinus floor elevation, typically 9-12 months should be allowed after implantation of bone graft material before placement of the titanium fixtures. X-rays should be taken to confirm the bone integrity prior to dental implant placement.

#### Periodontology

The filling of periodontal defects with *RegenerOss® Resorbable Xenograft* requires (along with plaque control) the successful local treatment of the periodontal lesion (e.g. root planing, debridement of granular tissue) prior to implantation.

#### Caution:

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

#### Adverse Reactions:

Possible complications that may occur with any dental surgery include swelling at the surgical site, flap sloughing, bleeding, local inflammation, bone loss, infection or pain.

#### How Supplied:

*RegenerOss® Resorbable Xenograft* is supplied sterile, non-pyrogenic, and for single use only.

Catalogue Number	Particle Sizes	Volume
ROXR05	0.25 - 1.0 mm	0.5 cc
ROXR10		1.0 cc
ROXR20		2.0 cc
ROXR40		4.0 cc
ROXRLG10	1.0 - 2.0 mm	1.0 cc
ROXRLG20		2.0 cc

#### Storage:

Store at room temperature (15°C/59°F - 30°C/86°F).

*RegenerOss®* is a registered trademark of Zimmer Biomet or its affiliates.

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