NOVABONE[®]
Dental Putty
Cartridge Dispenser

(US)

NOVABONE[®]
Dental Putty
Cartridge Dispenser

Instructions For Use

Indications:

The NovaBone Dental Putty Cartridge Dispenser is indicated for use only with NovaBone Dental Putty Cartridges. The cartridges are sold separately and come pre-filled with bone graft substitute.

(1) (2)

(3a) Correct (3b) Incorrect (4)

Instructions for use:

- Remove cartridge from the tray.
- Press the cartridge into the dispenser until it locks in place.
- Align the black plastic plunger and spring against the flat support on the underside of the handle lever, indicated by the arrow printed on the handle.
NOTE: Do not align with the short post on the handle as shown in Figure 3b, since the handle will not close or function properly in this position.
- Remove the small black cap from the tip of the cartridge.
- Squeeze the handle with a slow steady pressure to express the bone graft putty.
- Remove and dispose the cartridge after use.

Sterility:

The cartridge dispenser is provided non-sterile and will need sterilization prior to use. The dispenser can be heat sterilized either in an autoclave or using dry heat not exceeding 285° F (140° C), or by cold solution sterilization. Cleaning and sterilization instructions are provided below.

Contraindications / Precautions / Warnings:

See NovaBone Dental Putty directions for use as supplied with the device cartridges.

Caution:

Federal law (USA) restricts this device to sale by or on the order of a physician.

Cleaning and Sterilization Instructions

Cleaning:

At chair side after use, wipe of excess material and debris from cartridge dispenser.

Method 1: Disassemble unit by lifting rear handle up to allow removal of black plunger. Remove metal spring from plunger and place all three parts (handle, plunger and spring) into an ultrasonic cleaning unit. Cover and clean for 16-20 minutes (at 122°F/50°C if unit contains a heater). Remove and rinse thoroughly under sterile fresh warm water (de-mineralized water may be used). Shake to remove excess water and stand vertically to air dry.

Method 2: Disassemble unit as listed above. Thoroughly clean each part by hand scrubbing with a mild detergent and a soft bristle brush. Remove and rinse thoroughly under sterile fresh warm water (de-mineralized water may be used). Shake to remove excess water and stand vertically to air dry.

Sterilization:

The cartridge dispenser should be dry and fully reassembled before sterilization. Reassemble spring onto plunger, making sure the small spring hole is at back of plunger for proper reassembly (see Photo 1 to show correct method).

Photo 1

Next, insert the plunger and spring unit into the unit until the back of the plunger is in front of the platform on the underside of the rear handle lever (see Photo 2 to show correct method).

Photo 2

With the plunger in the correct position, the handle should move freely and not bind during use. If the plunger is not in the correct position, the handle will require excessive force to close, which may result in fracture of the handle or the plunger.

The cartridge dispenser may be sterilized either by autoclave or by cold solution sterilization.

Autoclave Sterilization: Proceed with the sterilization procedure according to the sterilizer manufacturer's instruction; normally, 20 minutes at 250°F/121°C – 15 psi.

Cold Solution Sterilization: Immerse the unit completely in activated Cold Solution at 25°C in a covered container for the length of time as recommended by the Cold Solution manufacturer. Remove unit from the solution and rinse thoroughly under sterile warm water (de-mineralized water may be used). Shake to remove excess water and stand vertically to air dry.

Manufacturer: NovaBone Products, LLC
 13510 NW US Highway 441
 Alachua, Florida 32615
 Tel: (386) 462-7660
 Fax: (386) 418-1636
 www.novabone.com

EU Authorized Representative: Peter Smith
 18 Yeates Close
 Thame, Oxfordshire, OX9 3AR UK

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