Gebrauchsinformation

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PUROS® ALLOGRAFT

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None

Form

Bone graft of human origin (allograft) as a pre-shaped graft or particles

Mode of action

For implantation. PUROS® ALLOGRAFT works as a biological scaffold for the new formation of resident vascularised bone. The process of resorption and new tissue formation generally begins after a few days post surgery and extends over several weeks and months. This remodeling process depends on the size of the bone graft and the remodeling capacity of the recipient site.

Indications

- Regeneration of periodontal bone defects
- Regeneration of furcation defects
- Regeneration after cyst resection or apicoectomy
- Regeneration of extraction sockets
- Regeneration of space between alveolar walls and dental implants
- Regeneration of defects after block harvesting
- Regeneration of gaps around block grafts
- Horizontal alveolar crest augmentation (particles)
- Sinus augmentation
- 3-dimensional (horizontal and/or vertical) alveolar crest augmentation (block augmentation)

Contraindications

- Active or latent infection in or around the implantation site.
- Any condition or diseases carrying an unacceptable increase in the postoperative risk.

Instructions for use

- Before implantation in patients, PUROS® ALLOGRAFT must be rehydrated with sterile 0.9% saline solution.
PUROS® ALLOGRAFT particulate material is covered by adding 0.9% saline solution. The graft material is stirred and allowed to stand for at least 2 min.
PUROS® ALLOGRAFT block grafts are placed into a sterile syringe with appropriate size. Sterile 0.9% saline solution is drawn up until the block graft is completely covered with fluid. The syringe is turned around and any air is expelled. The syringe nozzle is closed and the plunger is pulled to create a vacuum inside the syringe. The plunger is hold in this position to remove air out of the block and add saline solution. This process is repeated until no more air bubbles pass off the block (approx. 5-6 times).
Rehydrated PUROS® ALLOGRAFT may be stored in this solution throughout the surgical procedure until required. Rehydrated PUROS® ALLOGRAFT can be mixed with blood or blood components.
- The PUROS® ALLOGRAFT remains sterile as long as the sterile package is not opened and/or damaged. PUROS® ALLOGRAFT must be used before the expiration date. 
- Before clinical use the surgeon should thoroughly understand all aspects and limitations of surgical and grafting procedures. Patients should be informed about the limitations of the grafting procedures.

- In order to avoid dislocation and to allow stable bony incorporation of the graft material, PUROS® ALLOGRAFT should be applied to the surgical site in a stable manner. Ensure that the general principles for working in a sterile environment and for using PUROS® augmentation materials are followed.

Interactions

None known

Precautions

The Tutoplast® process is highly active against all kinds of pathogens and no cases of disease transmission have been confirmed in more than 1.500.000 surgical operations with Tutoplast® preserved tissue over the last 35 years. As with all biological products, the transmission of infectious diseases can not be ruled out categorically.

Adverse reactions

None known.

Patients should be instructed to inform the physician of any adverse reactions that may occur. Adverse reactions should be immediately reported to Zimmer Dental or Tutogen Medical GmbH.

Expiry date

PUROS® ALLOGRAFT has a shelf life of 5 years in the undamaged sterile package.


PUROS® ALLOGRAFT must not be used after the expiry date printed on the package and on the carton box. If the sterile package has been opened, but the PUROS® ALLOGRAFT has not been used, the bone graft must not be re-sterilized and must subsequently be discarded.

After opening the sterile package, any remaining PUROS® ALLOGRAFT must be discarded (not re-sterilisable).

Special precautions for storage

Store PUROS® ALLOGRAFT in a clean and dry place, protected from direct sun exposure at room temperature between 15°C and 30°C.

Special precautions for disposing of unused products

None

Date of the package insert

06/2017 "03"

Tutoplast® is a registered trademark of Tutogen Medical GmbH; PUROS® is a registered trademark of Zimmer Biomet or its affiliates.

Tutoplast® processed allografts have been successfully applied for more than 35 years in over 1.500.000 surgical procedures. This success is based, on one hand, on donor selection and serological tests. On the other hand, PUROS® ALLOGRAFT is preserved using the proprietary, multi-step Tutoplast® procedure that consists of meticulous cleaning of the tissue and gentle solvent dehydration monitored by continuous quality controls. This process removes cellular components, non-collagenous proteins and thus antigenic properties of the tissue. This process does not result in any chemical residues and is highly active against all types of pathogens. Preservation also includes sterilisation of the finished product by low-dose gamma irradiation. The natural composition of the bone matrix (collagen and calcium apatite) and consequently a high biomechanical stability of these allografts is preserved due to the Tutoplast® process.

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