RegenerOss®

Bone Graft Plug

INSTRUCTIONS FOR USE

A sterile, biocompatible bone graft composite plug for use in periodontal, oral and maxillofacial surgery.

DESCRIPTION

RegenerOss® Bone Graft Plug is composed of synthetic calcium phosphate granules and Type I collagen from bovine Achilles tendon in plug shaped forms designed for use in dental surgeries. The synthetic calcium phosphate mineral has a carbonate apatite structure, which is physically and chemically comparable to the mineralized matrix of human bone. The collagen component is highly purified type I collagen derived from bovine Achilles tendon. RegenerOss Bone Graft Plug is sterilized by gamma irradiation. The product is non-pyrogenic and for single use only.

PROPERTIES/ACTIONS

The mineral component of RegenerOss Bone Graft Plug has macro- and microscopic structures similar to human bone. The formation and ingrowth of new bone at the implantation site of RegenerOss Bone Graft Plug is favored, due to its trabecular architecture, interconnecting macro and micro pores and its natural consistency. The collagen facilitates handling of the bone mineral particles and acts to minimize migration of the granules.

The use of RegenerOss Bone Graft Plug may be considered when autogenous bone is not indicated, or insufficient in quantity to fulfill the needs of the proposed surgical procedure.

INDICATIONS AND USAGE

RegenerOss Bone Graft Plug is intended for use in dental surgery for:

- Augmentation or reconstructive treatment of the alveolar ridge.
- Filling of infrabony periodontal defects.
- Filling of defects after root resection, apicoectomy, and cystectomy.
- Filling of extraction sockets to enhance preservation of the alveolar ridge.
- Elevation of the maxillary sinus floor.
- Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR).
- Filling of peri-implant defects in conjunction with products intended for Guided Bone Regeneration (GBR).

CONTRAINDICATIONS

Contraindications customary to the use of bone grafts should be observed. RegenerOss Bone Graft Plug should not be used in patients with:


- Acute or chronic infection (osteomyelitis) at the surgical site
- Metabolic diseases (diabetes, hyperparathyroidism, osteomalacia)
- Severe renal dysfunction, severe liver disease
- High dose corticosteroid therapy
- Vascular impairment at the implant site

- Osteoporosis
- Known allergy or hypersensitivity to bovine derived implantable materials

ASEPTIC OPENING/PRESENTATION:

RegenerOss Bone Graft Plug products are packaged in a single layer blister tray. To ensure sterility is maintained, do not open the blister package until ready for implantation. Caution should be taken during opening, removal and implantation to make sure that the product does not come in contact with any non-sterile surfaces.

- To open:
 - Hold the blister tray in one hand with the top lid opening (square corner) away from you.
 - Pull the top lid back towards you.

- To remove:
 - Hold the flap away from both the contents of the package and the sterile field.
 - Drop the product onto the sterile field or if assisting with the procedure hold the exposed sterile item toward the clinician so he/she can grasp with a sterile instrument (e.g., forceps) or sterile gloves.
- To implant:
 - Prepare the device for implantation in accordance with the instructions for use below. Transfer the device from the sterile field to the patient defect area using a sterile instrument (e.g., forceps).

INSTRUCTIONS FOR USE

- After exposure of the bony defect with a mucoperiosteal flap, all granulation tissue must be carefully removed.
- If pre-hydration is desired, hydrate with patient's blood or sterile normal saline.
- In order to assure the formation of new bone, RegenerOss Bone Graft Plug should only be placed in direct contact with well vascularized bone. Cortical bone should be mechanically perforated.
- Loosely pack RegenerOss Bone Graft Plug into the osseous defect using a sterile instrument. Do not use excessive force, as this will result in compression of the particles and loss of trabecular architecture.
- Overfilling of the defects should be avoided.
- The mucoperiosteal flaps should be sutured to achieve primary closure, if possible. A surgical dressing may be placed over the wound for one to two weeks.
- If primary wound closure cannot be achieved completely, further immobilization of the flap (e.g., by incision through the periosteum) should be performed and/or a bioabsorbable membrane should be placed over the bone graft site.

PRECAUTIONS

In order to facilitate the formation of new bone, RegenerOss Bone Graft Plug should only be implanted in direct contact with a well vascularized bony tissue. Drilling may be recommended to facilitate bleeding from cortical bone. In larger defects a mixture of autogenous bone or bone marrow may improve the formation of new bone.

Implantology

Generally, in augmented areas, the placement of titanium fixtures should take place once the bone has sufficient strength and integrity for dental implant placement, which is typically greater than 6 months after implantation of a bone graft material. For sinus floor elevation, typically 9-12 months should be allowed after implantation of bone graft material before placement of the titanium fixtures. X-rays should be taken to confirm the bone integrity prior to dental implant placement.

Periodontology

The filling of periodontal defects with RegenerOss Bone Graft Plug requires (along with plaque control) the successful local treatment of the periodontal lesion (e.g. root planning, debridement of granular tissue) prior to implantation.

WARNINGS

RegenerOss Bone Graft Plug is not intended for load bearing applications. The plug shape is intended for ease of handling.

RegenerOss Bone Graft Plug cannot be re-sterilized or re-used. Open, unused RegenerOss Bone Graft Plug must be discarded. *In vivo* stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

Do not use if the product sterilization barrier or its packaging is compromised.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist.

ADVERSE REACTIONS

No adverse reactions have been reported in relationship to RegenerOss Bone Graft Plug. Due to the additional collagen component, allergic reactions cannot be totally excluded.

HOW SUPPLIED

RegenerOss Bone Graft Plug is supplied sterile, non-pyrogenic, and for single use only.

STORAGE

The product should be stored at room temperature (15°C/59°F to 30°C/86°F). Avoid excessive heat and humidity.

RegenerOss is a registered trademark of ZimVie or its affiliates.

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Manufacturer: Collagen Matrix, Inc. 15 Thornton Road Oakland, NJ 07436 USA
Distributed by: ZimVie Dental, LLC 4555 Riverside Drive Palm Beach Gardens, FL 33410 USA