## Dental Zirconia Blank Instruction for use

### 2. Software Design Notes

- Holding bars must be attached to each single-tooth restoration.
- In multi-unit restorations, end units must have 2 holding bars attached (oral and vestibular). Attach holding bars to other units as required.
- The diameter of the holding bars must be at least 2.0 mm.
- The holding bars must be attached at least 1.0 mm above the preparation edge.
- The holding bars should be placed at the anatomical equator of the tooth so that no under cuts are created and the restoration can be easily processed from both sides.
- In the case of multi-unit restorations with a pronounced curvature, a sintering support structure is advisable.
- Provide end units with a holding bar vertical to the sintering support structure.

#### 3. Milling

- The disc information can be transferred to any compatible CAM software.
- For manual entry the specific shrinkage factor can be found on the disc.
- When securing the disc in the holder, care must be taken to ensure that the circumferential notch and the disc holder are absolutely clean and that the screws are tightened evenly crosswise.
- Mill the blanks with the milling machine applying the milling strategies provided for the material.

- When milling the blanks, it is recommended to use cutters with sufficient cutting edge quality.

#### 4. Separating and finishing

- Suitable rotary tools (e.g. fine tungsten carbide burs) recommended for separating the restorations.
- Suitable rotary tools (e.g. fine tungsten carbide burs) are recommended to smooth out the attachment points of the holding bars. Rough tungsten carbide burs and/or grinding instruments are not suitable, as they may cause chipping, among other things.
- When using the brush infiltration technique, make sure the surfaces of the restorations are not contaminated with any grease or smudgy film since this may negatively influence the shade result.
- After finishing, thoroughly remove any zirconium oxide dust with a soft brush and then blast the restoration with oil-free compressed air.


#### 5. Optional: infiltration and drying

Infiltration technique is recommended recommended when using the unshaded blanks. The following notes should be taken into consideration:

- The coloring liquids must not be contaminated.
- The coloring liquids have to be sealed when not in use.
- If there is cloudiness or precipitation (e.g. sediments), the coloring liquids should no longer be used.
- Cloudiness or precipitation of the coloring liquids is promoted by contamination.
- Do not decant and/or store the coloring liquids in metal containers.
- Apply Coloring Liquids with a clean, metal-free brush on the restoration.
- Do not use Coloring Liquid Indicator in an unmixed state.
- Coloring Liquid Indicator colors are not stable for longer periods of time.
- Store mixed solutions in a sealed container and use them within 4 hours. After longer storage, the reproducibility of the tooth shade in the sintered state is no longer ensured.
- Infiltrated restorations must be completely dried before sintering.

#### 6. Sintering

Select a proper high-temperature sintering furnace to sinter the restorations. Please also refer to the operating instructions of the respective sintering furnace. The programming you use for your sintering furnace should be analogous to the sintering schedule below:

Sintering temperatures are recommendations. If necessary, carry out a trial sintering cycle and adapt the sintering temperatures or times as needed.

#### 7. Manual finishing after sintering

Sintered zirconia should be finished using diamond instruments under proper irrigation only. Keep the pressure on the framework material to a minimum and work only in one direction.

- Eliminate premature contacts one by one until the framework has reached its final position on the die(s).
- During the try-on and fitting of the framework, keep the dies on the cast and try on the framework as a whole.
- We recommend using a rubber polisher to smooth the basal side of the bridge connectors.
- Rinse the restoration under running water or use the steam jet to remove any adhering residue and glaze.

#### 8. Staining and drying

For individualization of full-contoured restorations, we recommend Upcera Glaze.

#### Processing at the dentists:

- Pre-processing: All prosthetic components must be cleaned, disinfected and/or sterilized prior to each use into the patient's mouth.
- Cementation: Dental Zirconia Blank restorations can be cemented conventionally with zinc phosphate cement or glass ionomer cement or adhesively using cement compounds. With conventional cementing, ensure sufficient retention and a minimum stump height of 4 mm. Please follow the operating instructions of the respective cementing.
- Post-processing: When grinding is required by the dentist while installing the restorations into the patient's mouth, the units must be well re-polished to a high gloss. For grinding, we recommend using only diamond grinding instruments. For polishing, we recommend using diamond polishing bits.

**Shelf period:** 5-year is an estimated replacement time. The product is not expected to degrade during this time period.

#### Storage:

- Keep container dry.
  - Store at room temperature.
- NOTE:** The manufacturer is not liable for any erroneous results since we have no control if this instruction is properly followed or if the product is appropriately used. Any claims for damages therefore shall be limited exclusively to the commercial value of our products.
- Please make sure that the dentist has received all the necessary information for the correct use of the product.
- Medical Device Adverse Event:** So far there is no medical device adverse event on this product.

**⚠ Cautions:** Federal Law restricts this device to sale by or on the order of a dentist.

#### Type and Class:

According to EN ISO 6872, the product belongs to Type II, Class4 or Class5. **SSCP:** According to the requirements of the Medical Device Regulation (MDR), SSCP and Basic UDI-DI of the products can be available in the European Medical Devices Database (Eudamed), website address: <https://ec.europa.eu/tools/eudamed>.

#### Safety instruction:

- Avoid collision with hard material, avoid squeeze or shaking.
- Don't apply to patient before the material is fully sintered.
- This product should only be handled by dental specialist.
- Wear suitable face mask to prevent fine dust inhalation. Avoid fine dust air pollution by using the air blaster in short bursts only rather than a continuous flow.

**Disposal:** The material is not hazardous to the environment.

#### European representative:

Name: MedNet EC-REP GmbH  
Address: Borkstrasse 10,48163 Münster, Germany.  
Distributor: Please refer to labeling.

#### Manufacturer:

Name: Liaoning Upcera Co., Ltd.  
Address: No.122 Xianghui Road, Economic Development Zone, Benxi, 117004 Liaoning, P.R. China  
Phone: +86 24-45565055  
Fax: +86 24-45854457  
E-mail: [ln.upcera@upcera.com.cn](mailto:ln.upcera@upcera.com.cn)

CE marking	Date of manufacture	Use-by date	Batch code	Non-sterile	Do not re-use
Model number	Unique device identifier	Medical device	Wear eye protection	Wear a mask	
Importer	Quantity	Serial representative			

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Naam: MedNet EC-REP GmbH
Naam: Liaoning Upcera Co., Ltd.