HT (High Translucent) & ML (Multi Layered) / HTML (High Translucent Multi Layered)

EN	ENGLISH	English	Dental Zirconia
FR	FRANÇAIS	French	Zircone dentaire
ES	ESPAÑOL	Spanish	Óxido de zirconio dental
IT	ITALIANO	Italian	Ossido di zirconio dentale
NL	NEDERLANDS	Dutch	Dentaaal zirkoniumoxide
DE	DEUTSCH	German	Dentales Zirkonoxid
SV	SVENSKA	Swedish	Dental zirkoniumdioxid
NO	NORSK	Norwegian	Dental zirkoniumdioksid
FI	SUOMI	Finnish	Hammaslääketieteellinen zirkoniumdioksidi
DA	DANSK	Danish	Dentalt zirkonia
PT	PORTUGUÊS	Portuguese	Óxido de zircónio dentário
EL	ΕΛΛΗΝΙΚΑ	Greek	Οδοντιατρικό οξειδίο του ζirkονίου
TR	TÜRKÇE	Turkish	Dental zirkonyum oksit
PL	POLSKI	Polish	Tlenek cyrkonu do aplikacji dentystycznych
RO	ROMÂNĂ	Romanian	Oxid de zirconiu dentar
HR	HRVATSKI	Croatian	Dentalni cirkonij
HU	MAGYAR	Hungarian	Fogászati cirkonoxid
SL	SLOVENŠČINA	Slovenian	Dentalni cirkonij
CS	ČESKY	Czech	Dentální zirkonoxid
SK	SLOVENSKY	Slovak	Dentálny oxid zirkónu
BG	БЪЛГАРСКИ	Bulgarian	Стоматологичен цирконий
UK	УКРАЇНСЬКА	Ukrainian	Стоматологічний цирконій
ET	EESTI	Estonian	Dentaalne tsirkoonium
LV	LATVIEŠU	Latvian	Zobārstniecības cirkonijs
LT	LIETUVIŠKAI	Lithuanian	Odontologinis cirkonio oksidas

MD	EN	Medical Device	DA	Medicinsk udstyr	CS	Zdravotnický prostředek
	FR	Dispositif médical	PT	Dispositivo médico	SK	Zdravotnícka pomôcka
	ES	Producto sanitario	EL	Ιατροτεχνολογικό προϊόν	BG	Медицинско изделие
	IT	Dispositivo medico	TR	Tıbbi cihaz	UK	Виріб медичного призначення
	NL	Medisch hulpmiddel	PL	Wyrób medyczny	ET	Meditsiiniseade
	DE	Medizinprodukt	RO	Dispozitiv medical	LV	Medicīniska ierīce
	SV	Medicinteknisk produkt	HR	Medicinski proizvod	LT	Medicinos priemonė
	NO	Medisinsk utstyr	HU	Orvosteknikai eszköz		
	FI	Lääkinnällinen laite	SL	Medicinski pripomoček		

Kuraray Noritake Dental Inc.

300 Higashiyama, Miyoshi-cho, Miyoshi, Aichi 470-0293, Japan
TEL +81-561-32-8953 FAX +81-561-32-8976

EMERGO EUROPE

Prinsessegracht 20
2514 AP The Hague, The Netherlands

Kuraray Europe GmbH

Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany
Phone: +49 (0)69 305 35835 Fax: +49 (0)69 305 98 35835
URL: <https://www.kuraraynoritake.eu>

HT (High Translucent) & ML (Multi Layered) / HTML (High Translucent Multi Layered)

CE 1639

I. Introduction

This IFU is for KATANA Zirconia HT (High Translucent) & ML (Multi Layered) / HTML (High Translucent Multi Layered). KATANA Zirconia HT & ML / HTML are a pre-sintered zirconia disc 98.5 mm in diameter, which contains a plastic ring. This is designed for all milling systems using this generic-type disc. (Please refer to your milling system's technical instructions for correct machine operation.) KATANA Zirconia HT has 6 available thicknesses: (10mm, 14mm, 18mm, 22mm, 26mm and 30mm) and has 3 shade variations: (HT10, HT12 and HT13). KATANA Zirconia ML has 3 available thicknesses: (14mm, 18mm and 22mm) and has 6 shade variations: (A Light, A Dark, B Light, C Light, D Light and A White). KATANA Zirconia HTML has 3 available thicknesses: (14mm, 18mm and 22mm) and has 14 shade variations: (A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D2, D3, NW). KATANA Zirconia HT is higher in translucency and KATANA Zirconia ML / HTML consists of 4 gradated shade layers. KATANA Zirconia HT & ML / HTML are recommended for use in fabricating FCZ (Full Contour Zirconia) restorations or the frameworks.


II. Intended Use

KATANA Zirconia is used for the fabrication of the all-ceramic restorations (frameworks, FCZ crowns, FCZ bridges, inlays, onlays and veneers).

III. Sintering Program


Sintering Program 1

Temperature	Programming Rate	Holding Time
Room Temp. — 1500°C (2732°F)	10°C/min. (18°F/min.)	—
1500°C (2732°F)	—	2 hrs
1500°C (2732°F) — Room Temp.	-10°C/min. (-18°F/min.)	—

Sintering Program 2

Temperature	Programming Rate	Holding Time
Room Temp. — 1515°C (2759°F)	35°C/min. (63°F/min.)	—
1515°C (2759°F)	—	0.5 hrs
1515°C (2759°F) — Room Temp. (*)	-45°C/min. (-81°F/min.)	—

IV. Composition

ZrO₂, Y₂O₃ etc.

V. Type and Class (ISO6872:2015)

Type:II/ Class:5

VI. Physical Properties

Coefficient of Thermal Expansion (25-500°C (77-932°F)): 9.9×10⁻⁶/K

VII. Directions for Use

- Take the disc from the packaging and confirm that the disc does not have a crack or other damage.
- Place the disc into the milling machine; then begin the milling process following the milling systems technical instructions.
- After milling, remove the restorations from the disc with a diamond bur, etc.
- The cutting waste or dust, which is attached to the restorations, can be removed with a gentle air stream.
- Put the restorations into the refractory sagger tray and place them in the sintering furnace.
- Depending on the performance of the sintering furnace used, review the sintering schedule shown above (III. Sintering Program) before sintering the restorations.
- After sintering, adjust the restorations with a diamond bur as needed.
- Confirm that the restorations have no cracks.
- 1 FCZ:
 - Glaze baking: Create a high shine surface by polishing, especially on the contact areas, then apply the glaze on all surfaces in the usual manner.
 - Hand polishing: Create a high shine surface on the entire restoration by polishing without using the glaze. In the case of using HTML, when the restoration is to be finished without the glaze, select a brighter shade than the final shade.
- 2 Frameworks: Build-up the dental porcelain (CERABIEN ZR or CZR PRESS LF) on the frameworks following the manufacturer's technical instructions. Check the coefficient of thermal expansion of the porcelain in the manufacturer's technical instructions to confirm compatibility.

VIII. Remarks on Handling

Contraindications:

- If the patient is hypersensitive to zirconia or any other components, this product must not be used.

Warning:

- If the patient or the dental professional demonstrates a hypersensitivity reaction, such as rash, dermatitis etc., discontinue use of the product and seek medical attention immediately.

Caution:

- This product should not be used when malocclusion, clenching or bruxism conditions are present.
- When milling the disc or cutting, grinding and polishing the restorations, use an approved dust mask and vacuum with air filter to protect your lungs from inhaling dust.
- When milling the disc or cutting, grinding and polishing the restorations, use safety glasses to prevent dust from getting into your eyes. If the dust gets into your eyes, immediately rinse with copious amounts of water and consult a physician.
- Do not use for any purposes except for dental restoration. This product is for dental application only.
- Do not touch the items heated in the furnace with your bare hands.
- Do not take the restorations out of the sintering furnace during high temperature, as the quenching causes the breaking. However, if an elevating type furnace is used under sintering program 2, the restorations may be removed from the furnace at 800 °C (1472°F) or less. When removed from the furnace, the restorations must be put on a tray made of ceramic fiber (for example: Noritake Porcelain Mat) to cool slowly.
- There may be a large difference between the setting temperature in the program and the temperature in the real furnace when using the furnace for sintering program 2. Please use the furnace after checking with the manufacturer that the furnace and temperature listed in sintering program 2 is compatible.
- Sintering program 2 is only recommended for frameworks (up to 3 units), FCZ crowns, FCZ bridges (up to 3 units), inlays, onlays and veneers.**
- The margins should be prepared with a deep chamfer and rounded shoulders, with cutting edges and corners rounded to eliminate sharp preparation corners. The angle of the axial surface should be within the range of 5 to 15 degrees.
- When preparing teeth, avoid the following: deep shoulders, J-margins, knife edges, serrated margins, non-tapered abutments, undercuts, guide grooves, the formation of retentive holes, and sharp corners.
- Keep the following thickness of this product for fabricating prosthetics:

Location & indication	Wall thickness
Anterior crown or bridge	0.4 mm or more
Veneer	0.4 mm or more
Posterior crown or bridge	0.5 mm or more
Inlay or onlay	0.5 mm or more

- Use the following cross-sectional areas for connectors when fabricating bridges.

Location & indication		Connector cross section
Anterior	2- or 3-unit bridges	7 mm ² or more
	more than 3-unit bridges	9 mm ² or more
Posterior	2- or 3-unit bridges	9 mm ² or more
	more than 3-unit bridges	9 mm ² or more

- Keep the connected bridge pontics maximum 2 teeth in bridge constructions. When the pontic of 2 teeth continues, keep the connector cross section between the pontic at 12 mm².
- Keep the cantilever bridge to 1 pontic tooth and the connector cross section at 12 mm².
- Choose a shade color that is brighter than the intended color for a thick restoration, as it may look duller depending on the thickness of the restorations.
- When using a sintering furnace for the first time and changing a sintering condition, colors after sintering may vary. Sinter a small piece of Zirconia beforehand and confirm the color.

Caution in conjunction with usage:

- Do not use the disc if there is a crack noticed after removing it from the package.
- Dry milling is recommended. If wet milling is used, the translucency of the disc may be reduced.
- If you find a crack in the restorations, do not use.
- The above sintering recommendation is only a guideline; some adjustments may be required depending on each individual furnace.
- When milling the disc, use caution when approaching the milling bur with the plastic ring to prevent detachment of the disc. Mill the disc as if leaving a zirconia part contacting the internal side of the plastic ring at 2 mm.
- Cut and remove the plastic ring prior to sintering the whole disc.
- Dispose of this product as a medical waste to prevent infection.

Storage:

- Store in a cool and dry place. Keep away from direct sunlight.
- The product should be stored at 10-30°C (50-86°F).
- Do not remove the disc from its packaging during storage.
- The disc is fragile, and requires care when handling.
- The product must be stored in an appropriate place where only dental personnel have access.
- The product must be used by the expiration date indicated on the package.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious accident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

"KATANA" and "CERABIEN" are trademarks of NORITAKE CO., LIMITED.

HT (High Translucent) & ML (Multi Layered) / HTML (High Translucent Multi Layered)

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