EN: English

Prosthetic Fixtures for Dental Implants

INSTRUCTIONS FOR USE

CAUTION – PLEASE READ CAREFULLY

Products distributed by Zimmer Biomet should only be used by dental specialists with experience in maxillofacial implantology and other specialties, such as dental diagnosis, planning, dental surgery, or prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. Products marked for single use exhibit a "Do not re-use" symbol on the product label. Products bearing this symbol on the label must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and / or other deterioration of the patient's health such as patient tissue infection.

All components distributed by Zimmer Biomet must be dry fitted before use to check that the fit is correct. The clinician is responsible for the correct application of restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please contact the manufacturer or your local distributor.

During any intraoral use and manipulation all products distributed by Zimmer Biomet must be secured to prevent aspiration due to its shape and size.

Caution: U.S. Federal law restricts this device to be sold to, or on the order of, a licensed dentist or physician.

We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must consist of a screw check. If the screws are subject to unusual wear, the complete integrity of the implant abutment should be checked. New screws must be used for any adjustment, correction, or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty.

INDICATIONS- PLEASE READ CAREFULLY

Implant abutments are used for prosthetic restorations of dental implants or for assisting procedures in the dental laboratory. The OverdenSURE® Attachment System is designed for use with overdentures or partial dentures, retained in whole or in part by endosseous implants in the mandible or maxilla.

CONTRAINDICATIONS – PLEASE READ CAREFULLY

All materials used are biocompatible; however, some patients may be allergic or hypersensitive to any of the materials and their components (specified in table 1).

Implant Attachment:

Not appropriate where a totally rigid connection is required. Use of a single implant with divergence of greater than 20 degrees from vertical is not recommended.

Attachments for Multi-Unit Abutments:

Not appropriate where a totally rigid connection is required.

COMPATIBILITY INFORMATION

All components distributed by Zimmer Biomet are available for a variety of connections (specified in table 2). For compatibility with dental implants and implant analogs please follow our catalog or guides or contact your local distributor.

(*) Some of the compatibilities are not cleared for United States. Please contact Zimmer Biomet for more details.

WARNINGS – PLEASE READ CAREFULLY

Products distributed by Zimmer Biomet should only be used by dental specialists with experience in maxillofacial implantology and other specialties, such as dental diagnosis, planning, dental surgery, or prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. Products marked for single use exhibit a "Do not re-use" symbol on the product label. Products bearing this symbol on the label must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and/or other deterioration of the patient's health like patient tissue infection.

Products distributed by Zimmer Biomet have not been evaluated for safety and compatibility in the Magnetic Resonance environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of products distributed by Zimmer Biomet in the Magnetic Resonance environment is unknown. Scanning a patient who has this device may result in patient injury.

PRECAUTIONS – PLEASE READ CAREFULLY

All components distributed by Zimmer Biomet must be dry fitted before use to check the correct fitting. The clinician will be responsible for correct application of those restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please call the manufacturer.

We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must include a screw check. If the screws have suffered unusual wear, the complete integrity of the implant-abutment should be checked. New screws must be used for any revision, correction, or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty. During any intraoral use and manipulation all products distributed by Zimmer Biomet must be secured to prevent aspiration due to their small size and shape

STERILIZATION AND RE-USE

All products distributed by Zimmer Biomet are supplied NON-STERILE. See table 1. The recommended sterilization cycle is standard gravity autoclave, exposure at 121 °C / 250 °F for 30 minutes with a drying time of 30 minutes (in accordance with the UNE-EN ISO 17665-1:2007 standard), using a sterilization wrap that is EU and FDA cleared for the indicated cycle.

Some devices are marked for "Single use only" because it is difficult or impossible to clean and decontaminate a used device, reuse can lead to cross-infection. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use device will not be accepted.

SHELF LIFE

The abutments are intended to be handled and manipulated by professional dental technicians before being shaped into their final form. They are labelled "for single use", but do not degrade over time and are not susceptible to variations in environmental or storage conditions. Therefore, neither the shelf life of the device nor the shelf life of the package is limited and no "use by" date is provided on the label.

Clinical screws are intended to be handled by professionals for implant-prosthesis fixation, and once assembled, they are not susceptible to variations in environmental conditions. Their shelf life is not limited. Therefore, a date is not labelled.

Table 1 General Information:

FAMILY	COMMERCIAL CODES	RAW MATERIAL	STERILIZATION	RE-USE
OverdenSURE® Abutment/ Class IIb	ODS-CERTxxx, ODS-EXHExxx, and ODS-TSVxxx	Titanium alloy ELI Ti-6Al-4V, ZrN coated	Autoclave before use on the patient	Single use
OverdenSURE® Abutment Retention Kit/ Class I	ODS-RIKIT2001 and ODS-RIKIT2002, ODS-RIKIT4001 and ODS-RIKIT4002, ODS-RIR4004, ODS-RIB2004, ODS-RIP2004, ODS-RIC2004, ODS-RIV4004, ODS-RIG4004, ODS-RH01, ODS-RH04, ODS-RH010, ODS-BS020	Retention Insert - Polynil (Nylon) Housing – Titanium Alloy (Ti-6Al-4V ELI) Blockout Spacer – White Silicone	The nylon Retention Inserts may be sterilized/ disinfected using an FDA approved liquid chemical sterilant. Nylon Retention Inserts must be soaked for a minimum of 3 hours in the liquid sterilant at room temperature	Single use
OverdenSURE® Analog	ODS-AA01, ODS-AA04, ODS-AA010	Surgical AISI-303	Autoclave	Not recommended
OverdenSURE® Coping	ODS-IC01, ODS-IC04, ODS-IC010	Surgical AISI-303 Screw: Titanium Alloy Ti-6Al-4V ELI	Autoclave	
OVERDENSURE® Bar Attachment	ODS-BARATT01, ODS-BARATT02, ODS-BARATT04	Titanium alloy ELI Ti-6Al-4V, ZrN coated	Autoclave before use on the patient	Single use
TOOLS / Class I	ODS-DRVR, ODS-IRTOOL	Handles: surgical stainless steel SUS 316L Tips: stainless steel AISI 420 MOD	Autoclave	In case of intraoral use: after sterilization

Table 2 Compatibility table.

COMPATIBLE WITH	PLATFORM	GINGIVAL HEIGHT
Zimmer Tapered Screw-Vent®	NP Green 3.5	1, 2, 3, 4, 5, 6mm
	NP Purple 4.5	1, 2, 3, 4, 5, 6mm
	WP Yellow 5.7	1, 2, 3, 4, 5, 6mm
3i External Hex	NP Purple 3.4	1, 2, 3, 4, 5, 6mm
	RP Blue 4.1	1, 2, 3, 4, 5, 6mm
	WP Yellow 5.0	1, 2, 3, 4, 5, 6mm
3i Certain®	NP Purple 3.4	1, 2, 3, 4, 5, 6mm
	RP Blue 4.1	1, 2, 3, 4, 5, 6mm
	WP Yellow 5.0	1, 2, 3, 4, 5, 6mm

APPLICATIONS AND DIRECTIONS FOR USE- OverdenSURE® Abutment

OverdenSURE® Abutment is a resilient implant attachment used for prosthetic restorations prepared by dental technicians in dental laboratories or by dental professionals chairside. Recommended for removable dentures.

Procedure: Measure the tissue thickness from the implant platform to the highest side of the gingival crest near the implant. Choose the corresponding height or 1.5 mm higher (the attachment should not be submerged) and put in place using the recommended torque from the Clinical Screw Recommended Torque document. The retention male part can be mounted in the dental clinic as well as in dental laboratory with self-curing resin using standard techniques.

APPLICATIONS AND DIRECTIONS FOR USE- OverdenSURE® Abutment Retention Kit

OverdenSURE® Abutment retention system is a part of the anchoring system for full or partial overdentures on implants, in cases where conventional full dentures have problems with retention and stability. The male retention kit can be mounted in the dental clinic using the self-curing resin, and in the dental laboratory, using the standard procedure.

Direct method: Place the silicone spacer over the attachment covering the entire area around it. If there is any space between the housing and soft tissue, put more spacers until the space disappears. Place the housing with the black Nylon male part on the OverdenSURE® Abutment Attachment in the patient’s mouth. Prepare a site for the housing in the prosthesis and perforate the buccal side to evacuate excess material. Apply the dual- or cold-curing resin around the housing and within the holes. Place the prosthesis in place while keeping the height of the occlusion and hold in position until the resin hardens. Remove the denture, spacers and cut off the excess. Polish before changing to the final male Nylon. **Indirect method:** same as direct, but on the master model with 14.054 analogues.

Retention force/divergence between implants: blue (1.5 lbs/680g, up to 20°), pink (3 lbs/1360g, up to 20°), clear (5 lbs/2270g, up to 20°), black (Processing Nylon male, 0 lbs, 0g, up to 40°), red (1 lbs/450g, up to 40°), orange (2 lbs/910g, up to 40°), green (4 lbs/1810g, up to 40°). ODS-RIKIT2001 and ODS-RIKIT2002 are kits up to 20°, ODS-RIKIT4001 and ODS-RIKIT4002 are kits up to 40°.

APPLICATIONS AND DIRECTIONS FOR USE- TOOLS

Tools are used for placement and adjustment of prosthetic restorations by dentist at the clinic or by dental technician in dental laboratory. Exchangeable screwdriver tips may be used for various implant systems. For intraoral use, all tools must be thoroughly autoclaved before use. Where using removable heads, make sure that it cannot come loose from the screwdriver. The screwdriver should always be fastened with a safety pin and cord attached to the groove designed for this purpose. It should be handled with care, taking the necessary precautions not to hurt the patient. Before using the screwdriver, check that it is the correct size and shape for the screw head that will be used. Only use screwdrivers that are in perfect conditions. A worn driver can make subsequent removal of the prosthesis impossible and can damage a screw head, so they should be regularly replaced.

In case of the need for adaptors, please assure that they correspond to 4 x 4 connection and to its abutment system (Multi-unit®, UniAbutment®, octagonal or Locator®).

SYMBOLS GLOSSARY

ANSI/AAMI/ ISO 15223-1:2012 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements.

Symbol	Title of Symbol (Reference Number)	
	Consult Instructions for Use	Indicates the need for the user to consult the instructions for use.
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Manufacturer	Indicates the medical device manufacturer.
	NON Sterile product	
	Medical Device	
	Medical product compliance marking controlled by the notified body number 0051	
	United States: Use only under prescription	

All products distributed by Zimmer Biomet are manufactured in accordance with ISO 9001, ISO 13485 and MDSAP requirements and bear the CE marking.

3i and 3i Certain® are registered trademarks of 3i Implant Innovations, USA

Zimmer Tapered Screw-Vent® is a registered trademark of Zimmer Dental INC.

Multi-unit® is a registered trademark of Nobel Biocare Services AG

UniAbutment® is a registered trademark of Dentsply Sirona

Locator® is a registered trademark of Zest Dental Solutions

These instructions replace all previous editions.

For any information about Products distributed by Zimmer Biomet, please contact your local distributor.

Distributed By:

Zimmer Biomet Dental

Global Headquarters

4555 Riverside Drive

Palm Beach Gardens, FL 33410

Phone: +1-561-776-6700

Fax: +1-561-776-1272

Zimmer Biomet Dental EMEA Headquarters

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Multi-unit® is een geregistreerd handelsmerk van Nobel Biocare Services AG