Universal Multi-Unit Prosthetic Fixtures for Dental Implants

INSTRUCTIONS FOR USE

PRODUCTS COVERED UNDER THIS IFU:

ZimVie Universal Multi-Unit Abutments, Abutment-level Ti-Bases and Components.

ELABELING:

The instructions for use can be accessed online by visiting the website provided on the label. Additional translations are also available in electronic format for download.

CAUTION – PLEASE READ CAREFULLY

Products covered in this document must be dry fitted before use to check that the fit is correct. The clinician is responsible for the correct application of restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt, please contact the manufacturer or your local distributor.

During any intraoral use and manipulation all products must be secured to prevent aspiration due to its shape and size.

Caution: U.S. Federal law restricts this device to be sold to, or on the order of, a licensed dentist or physician.

We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must consist of a screw check. If the screws are subject to unusual wear, the complete integrity of the implant abutment should be checked. New screws must be used for any adjustment, correction, or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty.

INDICATIONS- PLEASE READ CAREFULLY

Implant abutments are used for prosthetic restorations of dental implants or for assisting procedures in the dental laboratory.

The Universal Multi-Unit Abutment System is designed for multiple-unit screw-retained restorations with up to 30° of implant divergence where the following conditions apply:

- Adequate bone quality
- Adequate inter-arch space for a final hybrid framework
- Appropriate anterior/posterior spread no greater than 1.5mm for cantilever length
- Retained in whole or in part by endosseous implants in the mandible or maxilla.

Table 1 – Universal Multi-Unit Abutment Placement Indications

IMPLANT SYSTEM	PLATFORM	INTRAORAL POSITION	INDICATION
Eztetic™ Implant System, 3.1mm(D) TSX™ Implant System, 3.1mm(D)	NP -2.9mm(D)	Anterior	Yes
		Pre-molar	Splinted
		Molar	Not indicated
Tapered Screw-Vent® Implant System Trabecular Metal™ Implant System TSX™ Implant System	3.5mm(D) 4.5mm(D)	All	Yes
Tapered Screw-Vent® Implant System Trabecular Metal™ Implant System	5.7mm(D)	All	Yes

CONTRAINDICATIONS – PLEASE READ CAREFULLY

Do not use the titanium Ti-Base for restorations with cantilever on a single implant, bruxism, insufficient space, or direct metal-to-Ti-Base casting.

The 2.9mm platform Universal Multi-Unit Abutments should not be used in the molar region.

WARNINGS – PLEASE READ CAREFULLY

All materials used are biocompatible; however, some patients may be allergic or hypersensitive to any of the materials and their components.

Products distributed by ZimVie Dental, LLC should only be used by dental specialists with experience in maxillofacial implantology and other specialties, such as dental diagnosis, planning, dental surgery, or prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. Products marked for single use exhibit a "Do not re-use" symbol on the product label. Products bearing this symbol on the label must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and/or other deterioration of the patient's health like patient tissue infection. Do not re-sterilize product supplied sterile.

Products distributed by ZimVie Dental, LLC have not been evaluated for safety and compatibility in the Magnetic Resonance environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of products distributed by ZimVie Dental, LLC in the Magnetic Resonance environment is unknown. Scanning a patient who has this device may result in patient injury.

Magnetic Resonance (MR) Safety Information

MR Conditional

Warning: The RF safety of the device has not been tested. The patient may only be imaged by landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the RF coil.

A patient with this device can be scanned safely in an MR system under the following conditions:

Static Magnetic Field Strength (B0)	1.5 T and 3.0 T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	For body transmit coil, landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the coil. Extremity T/R coils permitted. Excludes Head T/R coil.
Operating Mode	Normal Operating Mode in the allowed imaging zone
Maximum Whole-Body specific absorption rate (SAR)	2 W/kg (Normal Operating Mode)
Maximum Head SAR	Not evaluated for head landmark
Scan Duration	No specific constraints due to implant heating

Note: removable restorations should be taken out prior to scanning.

STERILIZATION AND RE-USE

The Universal Multi-Unit Abutments covered in this document containing a "S" suffix are supplied STERILE. The multi-unit abutments are intended for single-use and should not be re-used nor reprocessed. They are subjected to a thorough manufacturing, control and cleaning process before being packaged in a sterile area and then further sterilized by gamma irradiation. Do not use them in case of package damage.

Some products covered in this document are supplied NON-STERILE. See Table 2. The recommended sterilization cycle is standard gravity autoclave, exposure at 121 °C / 250 °F for 30 minutes with a drying time of 30 minutes (in accordance with the UNE-EN ISO 17665-1 and 2), using a sterilization wrap that is EU and FDA cleared for the indicated cycle. Some Universal Multi-Unit Abutments covered in this document are supplied NON-STERILE. The non-sterile multi-unit abutments are intended for single-use and should not be re-used. They are subjected to a thorough manufacturing, control and cleaning process before being packaged.

Some devices are marked for "Single use only" because it is difficult or impossible to clean and decontaminate a used device, reuse can lead to cross-infection. Discard single-use and worn reusable instruments immediately after use. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use device will not be accepted.

Frequent use has a minimal effect on the instruments. The end of instrument life is usually determined by wear and tear and deterioration during use. Therefore, instruments can be reused with due care, provided they are not damaged or contaminated.

For the reuse of some instruments, they must be cleaned prior to sterilization. Here are some cleaning tips:

- Never place instruments of different types of materials together.
- Do not use metal brushes to remove impurities.
- Use disposable syringes for instrument cavities.
- When selecting detergents and disinfectants, make sure they are commercially available products and always follow the manufacturer's instructions.

SHELF LIFE

The abutments are intended to be handled and manipulated by professional dental technicians before being shaped into their final form. They are labelled "for single use", but do not degrade over time and are not susceptible to variations in environmental or storage conditions.

For the universal multi-unit abutments covered by this document that are supplied STERILE, read the labelling information carefully and do not use them beyond the expiry date informed on the label of their treatment (5 years from their date of manufacture).

Clinical screws are intended to be handled by professionals for implant-prosthesis fixation, and once assembled, they are not susceptible to variations in environmental conditions. Clinical screws and products supplied non-sterile do not have a shelf life. Therefore, neither the shelf life of the device nor the shelf life of the package is limited and no "use by" date is provided on the label.

Table 2 - General Information:

FAMILY	TYPE OF COMPONENT	REF CODE	RAW MATERIAL	STERILIZATION	RE-USE
MULTI UNIT	ABUTMENT	UMUA-NPXXXX UMUA-NPXXXX-17S UMUA-NPXXXX-30S UMUA-TSVXXXX UMUA-TSVXXXX-17S UMUA-TSVXXXX-30S	Titanium alloy ELI Ti-6Al-4V	Sterile - gamma	Single use
		UMUA-TSVXXX UMUA-TSVXXX-17 UMUA-TSVXXX-30	Titanium alloy ELI Ti-6Al-4V	Autoclave before intraoral use	Single use
MULTI UNIT COMPONENTS	HEALING CAP	UMUA-HCS	Titanium alloy ELI Ti-6Al-4V	Sterile - gamma	Single use
		UMUA-HC	Titanium alloy ELI Ti-6Al-4V	Autoclave before intraoral use	Single use
	IMPRESSION COPING	UMUA-CTICX UMUA-CTICXX UMUA-OTIC	Body: surgical stainless steel AISI-303 screw: Titanium alloy ELI Ti-6Al-4V		
	SCANBODY	UMUA-SBTB5, UMUA-SBIO	Polyether ether ketone (PEEK) Screw: Ti alloy ELI Ti-6Al-4V		
	DESKTOP SCANBODY	UMUA-SBDT	Polyether ether ketone (PEEK) Screw: Ti alloy ELI Ti-6Al-4V	No sterilization needed. Not indicated for intraoral use.	Reusable
	CYLINDER	UMUA-TC	Titanium alloy ELI Ti-6Al-4V	Autoclave before intraoral use	Single use
		UMUA-TCS	Titanium alloy ELI Ti-6Al-4V	Sterile - gamma	Single use
	ABUTMENT LEVEL TIBASE	UMUA-TTS UMUA-TTL UMUA-ASCT	Titanium alloy ELI Ti-6Al-4V	Autoclave before intraoral use	Single use
	SCREWS	UMUA-ASCTS UMUA-S	Titanium alloy ELI Ti-6Al-4V	Autoclave before intraoral use	Single use
	OVERDENSURE ATTACHMENT	UMUA-ODS	Titanium alloy ELI Ti-6Al-4V, ZrN coated	Autoclave before intraoral use	Single use
ANALOGS	UMUA-A UMUA-A10	Surgical stainless steel AISI-303	No sterilization needed. Not indicated for intraoral use.	Not recommended	
	UMUA-DA	Surgical stainless steel AISI-303	No sterilization needed. Not indicated for intraoral use.	Not recommended	
OTHER TOOLS	N/A	ASC-SD24 ASC-SD30 UMUA-D ISO4X4A	Handles: surgical stainless-steel SS 316L Tips: stainless steel AISI 420 MOD	Autoclave before intraoral use	Reusable

Table 3 - Compatibility Information.

UNIVERSAL MULTI UNIT ABUTMENTS

IMPLANT SYSTEM COMPATIBILITY	PLATFORM	ANGULATION	GINGIVAL HEIGHT
Eztetic™ Implant System, 3.1mm(D) TSX™ Implant System, 3.1mm(D)	NP - 2.9mm(D), Silver	Straight (0°)	1, 2, 3, 4, 5mm
		Angled (17° and 30°)	2.5mm, 3.5mm (17°) 3.5mm, 4.5mm (30°)
Tapered Screw-Vent® Implant System Trabecular Metal™ Implant System TSX™ Implant System	3.5mm(D), Green	Straight (0°)	1, 2, 3, 4, 5mm
		Angled (17° and 30°)	2.5mm, 3.5mm (17°) 3.5mm, 4.5mm (30°)
	4.5mm(D), Purple	Straight (0°)	1, 2, 3, 4, 5mm
		Angled (17° and 30°)	2.5mm, 3.5mm (17°) 3.5mm, 4.5mm (30°)
Tapered Screw-Vent® Implant System Trabecular Metal™ Implant System	5.7mm(D), Yellow	Straight (0°)	1, 2, 3, 4, 5mm
		Angled (17° and 30°)	2.5mm, 3.5mm (17°) 3.5mm, 4.5mm (30°)

APPLICATIONS AND DIRECTIONS FOR USE – ABUTMENT LEVEL TI-BASE

The abutment-level Ti-Base is used for prosthetic restorations prepared by dental technicians in a dental laboratory. The main use of the titanium Ti-Base is to support the bridge or individual zirconium dioxide elements made by CAD/CAM systems or with manual milling machines.

To use this product in a digital workflow, the ZimVie Universal Multi-Unit Abutment library and a scanbody are needed. To download and install one of our libraries please visit <https://www.zimvie.com/digital-libraries>.

Some dental CAD/CAM systems require a double scan with a diagnostic wax-up. If you need additional assistance with installation and/or proper use of the files and with use of the piece itself, please contact your local Territory Manager.

To improve cement adhesion, we recommend that the abutment level Ti-Base be thoroughly cleaned and degreased before cementation. You can use any implant luting cement (duo or chemical cured) extra-orally following the recommendations of the cement manufacturer. The internal surface of the zirconium structure (cementation zone) should be sandblasted and cleaned/degreased.

For secure grip, the diameter and height of the Ti-Base should not be reduced.

APPLICATIONS AND DIRECTIONS FOR USE – UNIVERSAL MULTI-UNIT ABUTMENTS

When placing straight Universal Multi-Unit Abutments, use plastic holder to facilitate insertion. When placing the angled Universal Multi-Unit Abutments use the metal holder to facilitate insertion in desired position.

Tighten the abutment using the Universal Multi-Unit Abutment tools in combination with the manual Torque Wrench (TWR), or the manual Low-Torque Indicating Ratchet Wrench (L-TIRW) and ISO Latch 4x4 Adapter (ISO4X4A) to a recommended torque of 30 Ncm.

Take an impression following the standard procedure (RevitaliZe Surgical and Restorative Technique Manual – ZBINST0023) and provisionalize the patient. If temporary prosthesis is not necessary, place the Universal Multi-Unit Abutment Healing Caps (UMUA-HC or UMUA-HCS).

APPLICATIONS AND DIRECTIONS FOR USE – HEALING CAP

Healing Caps are designed to cover the Universal Multi-Unit Abutment connection during the healing period.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSFER/ IMPRESSION COPINGS

The transfer or impression copings are used in closed and open tray techniques, for the transfer of the abutment position from the intraoral environment to the model in the dental laboratory. The transfer is used in the closed and open tray technique, for the transfer of the implant position from the intraoral situation to the model in the dental laboratory. Before being used, ensure that the connection seat for the implant is clean. Any traces of dirt could affect subsequent alignment of the prosthesis. Check compatibility with the implant model to which it will be connected. Remove the abutment, clean the connection with water/air and dry. Prepare the transfer in line with the implant connection and platform size. Place it on the implant and check that it is in the correct position. Put the screw in and tighten manually.

Before proceeding, check the final abutment seating using radiographic imaging.

Take an impression using the standard procedure (RevitaliZe Surgical and Restorative Technique Manual – ZBINST0023).

- The transfer is used in the closed and open tray technique, for the transfer of the implant position from the intraoral situation to the model in the dental laboratory. Before being used, ensure that the connection seat for the implant is clean. Any traces of dirt could affect subsequent alignment of the prosthesis. Check compatibility with the implant model to which it will be connected. Remove the abutment, clean the connection with water/air and dry. Prepare the transfer in line with the implant connection and platform size. Place it on the implant and check that it is in the correct position. Put the screw in and tighten manually. For the closed tray technique, after tightening lock the hex of the screw with wax. Use hydrocolloid, polyethylene or soft silicone. After hardening, place the transfer in its groove and check the stability before sending to the laboratory. This technique is strongly recommended only for implants without axial divergence.
- For the open tray technique choose the long screw. After tightening, position the impression tray such as to ensure full access to transfer screws from the outside. Before taking impressions, the transfers can be ferulised with acrylic resin placed on dental floss between consecutive implants. After hardening release all screws and remove the impression tray.

APPLICATIONS AND DIRECTIONS FOR USE – INTRA-ORAL AND DESKTOP SCAN BODIES

Desktop scan bodies are used for prosthetic restorations prepared by dental technicians in dental laboratories or by dentists at chairside in the CAD/CAM process. Designed to obtain geometric data from the master model using a desktop laboratory 3D scanner or for optical impressions using an intraoral 3D scanner.

For proper processing the Universal Multi-Unit digital library is needed. To download and install one of our libraries please visit <https://www.zimvie.com/digital-libraries>.

Before scanning the abutment and analog they should be visually checked for surface damage or implant connection for any tissue residue. For greater scanning precision, we recommend locating the flat surface of the scanning abutment in palatal/lingual orientation.

Fasten the abutment using the corresponding screw (directions in our catalogue and on the product label) by hand or with a maximum torque of 10 Ncm. In case of intraoral scan abutment, fix it with fingers using the embedded screw. The scanning abutment is a precision tool and being overtightened may change its geometry causing errors in the scanning process and discrepancy in accuracy. For most of scanners a scan-enhancing spray is not required.

The scanning process should be performed as recommended by the CAD/CAM system manufacturer. It is important to choose the correct connection in the software and the corresponding type for the chosen restoration (engaging/non-engaging). After the process, the scanbody can be loosened and placed gently on the tray or box. To complete the scan, gathering additional information (e.g., occlusion index, silicone bite, gingiva shape) may be required. In case of intraoral use, sterilization is required.

APPLICATIONS AND DIRECTIONS FOR USE – CYLINDER

Cylinders are used as support for provisional screwed crowns, bridges and prosthesis.

APPLICATIONS AND DIRECTIONS FOR USE – CLINICAL SCREW

The screw is for fixing prosthetic restorations and auxiliary abutments over the Multi-Unit Abutment or analog. Make sure to torque the Universal Multi-Unit Abutment Screw (UMUA-S) to 15 Ncm. For best results, the following conditions must be adhered to meticulously:

- Use the suitable model screwdriver and size for tightening and unscrewing. If in doubt, check that the next size screwdriver does not fit into the seat. The screwdriver should be placed on the longitudinal axis of the prosthesis/implant assembly.
- A new screw should be used when assembling a prosthesis for the first time and for every check thereafter.
- When transferring to the patient, use a new screw.
- **Position the patient to avoid aspiration should the screw fall during screwing/unscrewing.**

APPLICATIONS AND DIRECTIONS FOR USE – OVERDENSURE ATTACHMENT

OverdenSURE attachment is a resilient attachment used for use with the universal multi-unit abutment. The removable denture restoration is prepared by a dental technician in a dental laboratory or by a dentist chairside.

APPLICATIONS AND DIRECTIONS FOR USE – ANALOGS


The analogs are for simulating the multi-unit abutment connection and position on a stone master model. Before any laboratory manipulation ensure that the anti-rotation and anti-removal parts on the analog are securely affixed. Check that the analog and the prosthetic element connections match in size and type before tightening. A sealed and passive connection is needed. Do not re-use an analog, because such a connection may be impossible to achieve with a product that may no longer conform to specifications.

APPLICATIONS AND DIRECTIONS FOR USE- TOOLS

Tools are used for placement and adjustment of prosthetic restorations by dentist at the clinic or by dental technician in dental laboratory. For intraoral use, all tools must be thoroughly sterilized before use. Tools should be handled with care, taking the necessary precautions not to hurt the patient. Before using the screwdriver, check that it is the correct size and shape for the screw head that will be used. Only use screwdrivers that are in perfect condition. A worn driver can make subsequent removal of the prosthesis impossible and can damage a screw head, so they should be regularly replaced.

SYMBOLS GLOSSARY

ANSI/AAMI/ ISO 15223-1:2012 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements.

Symbol	Title of Symbol (Reference Number)	
	Consult Instructions for Use or Consult Electronic Instructions for Use ifu.biomet3i.com	Indicates the need for the user to consult the instructions for use.
	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Manufacturer	Indicates the medical device manufacturer.
	Date of Manufacture	Indicates the date when the medical device was manufactured for non-sterile products
	Non Sterile product	
	MR Safety	Magnetic Resonance (MR) Safety Information
	Medical Device (applicable only for non-sterile products)	
	Medical product compliance marking controlled by the notified body number 0051 except sterile Multi unit abutments (see label information for each case)	
Rx Only	United States: Use only under prescription	
	Sterilized using irradiation	
	Use by date	
	Do not use if package is damaged	
	Do not re-sterilize	

	<p>Indicates a single sterile barrier system</p>
	<p>Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences</p>

All products distributed by ZimVie Dental, LLC are manufactured in accordance with ISO 9001, ISO 13485 and MDSAP requirements and bear the CE marking (non sterile products only).

These instructions replace all previous editions.

For any information about Products distributed by ZimVie Dental, LLC, please contact your local Territory Manager.

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„Biomet 3i, LLC“
„Biomet 3i, LLC dba ZimVie Dental, LLC“