Instruction for use

Zfx™ Zirconium Color Effect

Version: 4/09.2021
Compiled at 26.07.2021
Effective date: 01.09.2021
Replaces version: 3/04.2018
Page 1 of 2

Material-related characteristics

Chemical Composition [wt%]

ZrO ₂ + HfO ₂ + Y ₂ O ₃	≥ 99
Al ₂ O ₃	≤ 0.1
Other oxides	< 0.1

Physical specifications

Density (after sintering)	[g/cm ³]	> 6.0
CTE (25°C – 500°C)	[10 ⁻⁶ K ⁻¹]	~ 10.5
Fracture toughness (SEVNB)	[MPa·m]	> 10
Fracture toughness (SEPB)	[MPa·m]	4.0
Flexural strength	[MPa]	1200 (± 250)
E modulus	[GPa]	> 200

Explanation of symbols

Manufacturer Information: Address of legal manufacturer facility

Date of manufacture

Distributor Information: Address of legal distributor facility

Use by

LOT-Number, Symbol followed by Production Date of the lot.

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

Store dry

Consult instructions for use
www.zfx-dental.com

CE Mark for Medical Devices Class 2

Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

Instruction Manual

1. Indications

Zfx™ Zirconium Color Effect are dental mill blanks made from zirconium dioxide for the manufacture of dental prostheses. Zfx™ Zirconium Color Effect can be used for all dental constructions incl. bridges >3 units with up to two connected pontics in the front and side tooth area.

2. Contraindications

- x Parafunctional habits
- x insufficient space
- x inadequate preparation
- x known intolerance to the contained components
- x insufficient oral hygiene

3. Safety information

Please pay attention to the information in the safety data sheet in its current version.

Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

4. Handling and storage

The blanks are manufactured and tested to highest quality standards. Verify the integrity of the packaging and the blank before the first processing. Check if the content of the packaging corresponds to the declaration given on the label. Do not use damaged material. Storage only in the original containers in a cool and dry environment. Avoid vibrations, contaminations and contact with fluids.

5. Processing / Designing

Handling this medical device should only be allowed to trained staff.

Generally, the following construction parameters need to be considered:

Indication	Scheme	Wall thickness [mm]	Connector cross section [mm ²]
Single crown	X	0.4	-
Telescope	X	0.5	-
Bridge 3 units Anterior / Posterior	z.B. XOX	0.5/0.5	>7/>9
Bridge 4 or more units Anterior / Posterior	z.B. XOOX XOOXOOX	0.6/0.6	>9/>12

Legend: X = abutment tooth O = pontic

Please note: Depending on the construction, the connector cross section might need to be of larger dimensions. An oval connector cross section is desirable; the height of the connector is the decisive factor. Frameworks for ceramic veneering need to be designed in a way to support the veneering ceramic in the cusp region and allow a constant layer thickness. Shoulder or chamfer preparation is recommended.

6. Milling, sintering and further processing

Machining of blanks should always be conducted in the corresponding milling system. All information given by the manufacturer of the machine shall be taken into account. Information for Amann Girrbach users: Please find code (scale factor) for utilization of 98 mm discs in Amann Girrbach milling machines on blank margin. The frameworks need to be examined for visible defects.

Instruction for use

Zfx™ Zirconium

Color Effect

Version: 4 / 09.2021
Compiled at 26.07.2021
Effective date: 01.09.2021
Replaces version: 3 / 04.2018
Page 2 of 2

Do not process damaged or contaminated frameworks. The white framework can be colored before final sintering with Zfx™ Color Liquids Allround – AoS (pay attention to separate instruction manual).

Please pay attention to our separate sintering instruction!

Sintering cycle for normal furnace filling without cover:

- ↑ Heating up to 900°C (1652°F) with 8°C / min (46°F / min),
- dwell at 900°C (1652°F) for 30 minutes,
- ↑ heating up to final temp. 1450°C (2642°F) with 3°C / min (37°F / min),
- dwell at final sintering temperature 1450°C (2642°F) for 120 min,
- ↓ cooling to 200°C (392°F) with 10°C / min (50°F / min).

Avoid additional mechanical actions after sintering like blasting or grinding. If corrections are inevitable use water-cooled tools for the conditioning and make sure that there is no development of heat, which could lead to cracks in the material. Please work only with sharp, diamond-coated tools and very slight contact pressure. Areas that are under tension in clinical use (e.g. connectors) must not be processed. Do not separate in interdental regions. Avoid sharp edges.

7. Ceramic veneering

Please use a veneering ceramic with a suitable WAK and pay attention to the manufacturer's recommendation. Slowing down heating and cooling rate in case of massive constructions is strongly recommended.

Weight per tooth unit [g]	< 1	2	3	> 4
Heating & cooling rate °C / min [°F / min]	55 [131]	45 [113]	35 [95]	25 [77]

8. Mounting

For mounting we recommend conventional cementation with zinc oxide phosphates cements or glasionomer cements. Also luting composites may be used. Sufficient retentions and a stump with height of at least 3 mm are essential. A temporary mounting is not recommended!

9. Information on Manufacturer

Dental Direkt GmbH
Industriezentrum 106 - 108
32139 Spenge
Germany
T +49 (0) 5225 / 86 319 - 0
F +49 (0) 5225 / 86 319 - 99
info@dentaldirekt.de
www.dentaldirekt.de

CE 0482

10. Information on Distributor

Zfx GmbH
Kopernikusstraße 15
85221 Dachau
Germany
T +49 (0) 8131 / 33 244 - 0
F +49 (0) 8131 / 33 244 - 10
info@zfx-dental.com
www.zfx-dental.com

We reserve the right to make changes as a result of the continuous development of our products. Please find the current version of the instruction manual on our website: www.zfx-dental.com

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