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
Instruction for use - English

High & Low Torque Indicating Ratchet Wrench

1 Intended purpose

The Low Torque Indication Ratchet Wrench is used for the adjustment of abutments, abutment screws and prosthetic screws during oral surgery and prosthetic procedures. The High Torque Indicating Ratchet Wrench is used for placement and adjustment of dental implants during oral surgery.

2 Product Description

The Torque Wrench (main body, head and Driver) is made of lightweight titanium and stainless steel. The Torque Wrench & Driver product design, finish and quality make the tool well suited for placement and adjustment of implants, abutments and screws during surgery and other clinical procedures. The Torque Wrench & Driver is easy to prepare for cleaning and sterilization. Exchangeable screwdrivers and the Driver make the tool compatible to various implant systems. The driver is available for a variety of connections. For specific product descriptions, please refer to individual product labels.

3 Chemical composition

Chemical composition of materials used for the Torque Wrench and Driver is presented in table 1-2.

Table 1: Chemical composition of the Torque Wrench. Content of below-mentioned materials depends on the variant.

Material	Composition	Patient contacting	CAS nr.
Titanium alloy (Ti-6Al-4V-ELI)	90 wt.% Ti, 6 wt.% Al and 4 wt.% V	Yes (Pure Metal)	99906-66-8
Sandvik 4C27A	84.8 wt.% Fe, 13 wt.% Cr, 1.2 wt.% Mo, 0.8 wt.% Ni, 0.2 wt.% S	No (internal component)	65997-19-5
AISI 304	67-71.5 wt.% Fe, 17.5-19.5 wt.% Cr, 8.0-10.5 wt.% Ni, 2 wt.% Mn, 1 wt.% Si	Yes (Pure Metal)	65997-19-5
Tungsten carbide and Amorphous carbon (Diamond like carbon)	80 at.% C + 16.8 at.% W with residual Ni	Yes (Outer layer)	12070-12-1 (WC) + 7782-42-5 (C)
Tungsten carbide binding layer	50 at.% W + 50 at.% C (WC)	No (Inner layer)	12070-12-1

Table 2: Chemical composition of the Driver. Content of below-mentioned materials depends on the variant.

Material	Composition	Patient contacting	CAS nr.
AISI 304	67-71.5 wt.% Fe, 17.5-19.5 wt.% Cr, 8.0-10.5 wt.% Ni, 2 wt.% Mn, 1 wt.% Si	Yes (Pure Metal)	65997-19-5
Sandvik 4C27A	84.8 wt.% Fe, 13 wt.% Cr, 1.2 wt.% Mo, 0.8 wt.% Ni, 0.2 wt.% S	Yes (Pure Metal)	65997-19-5
Sandvik 1RK91	71.7 wt.% Fe, 12 wt.% Cr, 9 wt.% Ni, 4 wt.% Mo, 2 wt.% Cu, 0.9 wt.% Ti and 0.4 wt.% Al	Yes (Pure Metal)	65997-19-5
Tungsten carbide and Amorphous carbon (Diamond like carbon)	80 at.% C + 16.8 at.% W with residual Ni	Yes (Outer layer)	12070-12-1 (WC) + 7782-42-5 (C)
Tungsten carbide binding layer	50 at.% W + 50 at.% C (WC)	No (Inner layer)	12070-12-1

4 Basic UDI

Torque Wrench 5712821001001U2
Driver 5712821002001U9

5 Indication for use

The interlinked Torque wrench and driver are adaptable to a variety of tools and together they apply a specified amount of torque to a tool. None of the tools are active devices.

The torque indicated on the torque wrench should correspond to that specified by the manufacturer of the implant, prosthetic screw or other component to be installed. The driver used should be compatible with the tool to be used to install the implant, prosthetic screw or ISO 1797 connection.

6 Intended patient population

The Torque Wrench & Driver is intended to be used on patients that are considered eligible for dental implant treatment. No special training of the patient is required.

7 Intended users

Primary users (during surgery and dental procedures): dentists, dental technicians, and surgeon professionals
Secondary users (cleaning, storing etc.): dental and surgeon assistants.

Use environment

Dental laboratories, dental clinics or operating room.

8 Clinical benefits

The Torque Wrench and Driver does not have a clinical benefit. The Torque Wrench and Driver can successfully assist as a tool to be used to install the implant, prosthetic screw or other component during dental surgeries.

9 Expected lifetime

Expected lifetime of up to 5 years.

10 Use

The Torque Wrench & Driver is delivered non-sterile. Prior to use it must be cleaned with water and brush followed by sterilization according to section 12. The use of the Torque Wrench & Driver is described below:

- After cleaning make sure the first line on the scale of the Torque Wrench arm aligns with the arrow (see figure 1.A).
- Attach the head to the body of the Torque Wrench by pushing the components together and turning them in opposite direction until an audible click is achieved (see figure 2).
- The arrow on the head of the Torque Wrench shows the direction in which the Torque Wrench & Driver is functioning, and the head is labeled "IN" (see figure 3.A) and "OUT" (see figure 3.B). "IN" denotes clockwise rotation; "OUT" denotes counterclockwise rotation. In order to change direction, flip the Torque Wrench to the opposite side (see figure 3) and push the Driver into the head of the Torque Wrench until an audible click is achieved (see figure 4).
- Insert the tool into the Driver.
- Place the tool in the implant/screw and turn the Torque Wrench arm in the direction of the arrow until the desired torque is achieved (see figure 6). The applied torque must be in accordance with the screw/implant manufacturers' instructions.
- After each use clean and sterilize the Torque Wrench & Driver according to section 12.
- Inspect instruments under normal lighting for the removal of visible soil and for damage or wear: Any instrument showing visible damage or wear must be disposed of according to section 16.

11 Warnings and precautions

Potential adverse occurrences may include but are not limited to the following list. Proper placement and use of this instrument may help to mitigate these risks.

- If overloaded or dropped**

If overloaded, dropped or in other way mishandled, the product must no longer be used since correct functionality can no longer be guaranteed.

- Do not use hydrogen peroxide**

This product must be sterilized prior to use and it must not be cleaned by using hydrogen peroxide.

- Scale alignment**

The Torque Wrench & Driver does not require recalibration but before each use, make sure that the Torque Wrench & Driver is fully functional and that the first line on the scale aligns with the arrow (see figure 1.A).

- Ensure correct assembly**

Before each use, make sure that all parts are correct assembled (no loose parts are allowed).

- Do not exceed max scale**

The arm of the Torque Wrench must not go beyond the end of the scale (see figure 6.B), as this could result in permanent deformation of the Torque Wrench arm i.e. future inaccurate torque readings (see figure 1.B).

- **No scale indication when used as ratchet wrench**

If the Torque Wrench & Driver is used as a ratchet wrench it does NOT provide the user any indication of what torque that is used.

Notice: When the Elos Torque Wrench & Driver is used as a ratchet wrench, specified max torque on Driver must not be exceeded. If no max torque on Driver is specified, the torque applied must not exceed 150Ncm.

- **Use screw/screwdriver/Driver/implant manufacturer torque recommendations**

Please check the individual screw/screwdriver/ Driver/implant manufacturers' instruction for use or label regarding the maximum allowed torque. Do not under any circumstances apply higher torque than stated even if the Torque Wrench & Driver is used as a ratchet wrench (without using the Torque Wrench arm).

- **Cleaning – do not use metal brush or steel wool**

Do not use a metal brush or steel wool to clean the Torque Wrench or Driver as this will damage the instrument.

- **Only for use by professionals**

The Torque Wrench & Driver must only be used by professionals.

12 Cleaning and sterilization

The following cleaning and sterilization method is recommended. Further instruction can be found in "Cleaning and sterilization guideline" at <https://elosmedtech.com/IFU/>

1. Disassemble the Torque Wrench (see figure 5)
2. The separated parts are now ready to be cleaned according to the "Cleaning and sterilization guideline".
3. Sterilization according to the "Cleaning and sterilization guideline". Summary below:

Note: For US: Steam Sterilization Cycle 132°C (270°F)/4min.

Procedure	Dynamic-air-removal sterilization Cycle	
Exposure time	4 min.	3 min.
Temperature	132°C (270°F)	134° C (273°F)
Drying time	20 min.	

Caution: Accompanying packaging material CANNOT be autoclaved unless it is explicitly stated that the material is suitable for sterilization.

13 Further information

For additional information about the use of the Torque Wrench & Driver, please contact your local sales representative.

14 Validity

Upon publication of this instruction for use, all previous versions are superseded.

15 Storage and Handling

Recommended storage temperature range is 10-40°C (50-104°F) and the recommended relative humidity range is 0-80%

16 Disposal

When disposed of the Torque Wrench & Driver may present an infection or microbial hazard and should therefore be disposed of as biological waste, according to national requirements.

17 Serious Incidents:

If any serious incidents occur in relation with this device, they should be reported to Elos Medtech Pinol A/S and to the competent authorities in the country where the incident occurred. A serious incident is a situation where the device directly or indirectly led or might have led to death, a serious deterioration in the health of a person or a serious public health threat. Contact Elos Medtech Pinol A/S at: complaint.empi@elosmedtech.com

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist or physician.

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complaint.empi@elosmedtech.com adresinden Elos Medtech Pinol A/S ile temasa geçin