




**Step 4:** Remove the ISRT07 Manual Reversing Drill and then remove the ISRT05 Guide Handle from the implant.



**Step 5:** Insert the ISRT06 Screw Removal Tool into the implant and seat the tool onto the top portion of the broken screw. Press onto the screw with slight to moderate pressure to engage/capture the screw. Rotate the tool in the reverse direction (counter-clockwise) several rotations. The screw should release from the implant internal threads.



**Step 6:** Pulling the ISRT06 Screw Removal Tool from the implant. The fractured screw should remain captured in the ISRT06. The screw can be removed by pulling it out of the tool tip. If the tool fails to catch the screw, repeat steps 2 – 5 until the screw is extracted from the implant.



**Step 7:** After the broken screw has been removed, insert the IWSU30 Certain Waxing Screw/Guide Pin to evaluate the integrity of the internal threads of the implant. If the IWSU30 does not rotate easily into the implant, significant thread damage may have occurred to the implant threads and the implant may need to be removed. If the IWSU30 rotates easily into the implant, you may proceed with placement of the new abutment.



Batch code



Caution, consult accompanying documents



Consult Instructions for Use  
[www.ifu.biomet3i.com](http://www.ifu.biomet3i.com)



Use by



Legal Manufacturer



Date of Manufacture



Catalogue number



By prescription only



Authorized European Representative



Do not use if package is damaged



Non-sterile



























