Instructions for Use

CopiOs Extend™ Membrane

The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.

CopiOs Extend Membrane should be essentially resorbed in approximately 6 to 9 months following surgery. However, probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

Distributed by:

Zimmer Biomet Dental
4555 Riverside Drive
Palm Beach Gardens
FL 33410 USA
+1 561.776.6700

zimmerbiometdental.com

ENGLISH

Instructions for Use

CopiOs Extend™ Membrane

Before using this product, the surgeon/practitioner should carefully study the indications, contraindications, recommendations, warnings and instructions and fully comply with them. The manufacturer, the importer and the distributor of these products are not liable for complications, other negative effects or damages that might occur for reasons such as incorrect indications or surgical technique, unsuitable choice of material or handling thereof, unsuitable use or handling of the instruments, asepis and so on. The surgeon/practitioner is responsible for any such complications or other consequences. It is also the surgeon/practitioner's responsibility to properly instruct and inform the patient on the functions, handling and necessary care of the product and on all known product risks.

Description

CopiOs Extend Membrane is a white, nonfriable membrane matrix engineered from highly purified porcine dermis. CopiOs Extend Membrane is resorbable in 6 to 9 months which eliminates the need for a second surgical procedure that is normally required to remove a non-resorbable membrane.

Indications

CopiOs Extend Membrane is intended for use in oral surgical procedures as a resorbable material for:

- Use in augmentation around implants placed in immediate and delayed extraction sockets
- Localised ridge augmentation for later implantation
- Alveolar ridge reconstruction for prosthetic treatment
- Filling of bone defects
- Guided bone regeneration in dehiscence defects
- Guided tissue regeneration procedures in periodontal defects

Contraindications

CopiOs Extend Membrane is contraindicated in patients who have:

- Acute infections or contaminated wounds in the oral cavity
- Known allergy to collagen of animal origin or to other porcine-derived products
- Clinically significant renal, hepatic, cardiac, endocrine, hematologic, autoimmune or systemic disease, which in the physician's judgment, will prevent safe implantation or likely healing

Warnings

Clinicians should use care in screening their patients for any known allergies to collagen or porcine-derived products.

Precautions

As with all surgical procedures, caution should be exercised when treating medically compromised patients such as patients receiving long-term steroidal therapy or currently taking anticoagulants. Patients with clinically significant systemic diseases, indicating a history of anaphylactic reactions, autoimmune diseases, uncontrolled diabetes or severe hypertension have not been implanted with the membrane; therefore, the safety and effectiveness for those patients have not been determined. The membrane has not been evaluated in pregnant women, children and/or patients with conditions involving extremely severe defects with little periodontium or bone.

CopiOs Extend Membrane cannot be re-sterilized or reused. Open, unused CopiOs Extend Membrane must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if reused.

Do not use if the product sterilization barrier or its packaging is compromised.

Adverse Reactions

Possible complications that can occur with any dental surgery include infection, swelling of the intraoral tissue, thermal sensitivity, gingival recession, excessive gingival bleeding, flap sloughing, resorption or ankylosis, with loss of crestal bone height, pain, or complications associated with the use of anesthesia. Minor discomfort may occur for a few days.

Sterility

The product is supplied sterile, non-pyrogenic and for single use only.

Administration

CopiOs Extend Membrane is packaged in a double sterile pouch. The outer pouch should be opened carefully, allowing the inner pouch to be placed onto a sterile field. The membrane should be removed from the inner pouch with sterile gloves or instruments.

The bone defect is exposed by a mucoperiosteal flap and basic surgical procedures are performed (e.g. curettage). The space-making material such as autologous bone, mineralized bone particles, demineralized bone matrix and inorganic grafting materials may be used to fill the defect.

CopiOs Extend Membrane can be placed either dry or hydrated. If the clinician prefers the handling characteristics of the hydrated collagen, the membrane can be hydrated in sterile water or saline solution for approximately five minutes prior to the final placement.

CopiOs Extend Membrane can be trimmed to the size and shape of the defect in the dry or wet state using sharp, sterile scissors.

CopiOs Extend Membrane should overlap the walls of the defect by at least 2mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material.

Fixation of the membrane may be indicated to avoid displacement due to loading or mobilization. The membrane can be sutured in place using absorbable sutures and a non-cutting needle. Resorbable tacks can also be used to affix the membrane. The mucoperiosteal flap is sutured over the collagen membrane and the wound should be closed completely to avoid accelerated resorption due to membrane exposure. The membrane is expected to be fully resorbed in 6 to 9 months.

Post-operative Procedures

CopiOs Extend Membrane is completely resorbable and should not be removed. Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for four weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent.

The patient should refrain from brushing the treated area for two weeks following the surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Dental floss should not be used prior to four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated.

The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.

CopiOs Extend Membrane should be essentially resorbed in approximately 6 to 9 months following surgery. However, probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

Storage

The product should be stored at room temperature. Avoid excessive heat and humidity.

How supplied

One (1) membrane per package, available in a variety of sizes:

15mm x 20mm; 20mm x 30mm; or 30mm x 40mm

Caution

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician or dentist.

©2017 Zimmer Biomet. All rights reserved.

CopiOs Extend is a registered trademark of Zimmer Biomet or its affiliates.

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